CysC (GFR) control kit

for quality control procedures in examining the accuracy and precision of the quantitative in vitro determination of Cystatin C with CysC (GFR) test kit on smart or CUBE laboratory photometer.

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Indication
CysC (GFR) control kit

Kit size
2 x 1 ml / low and high

Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

The CysC (GFR) control kit is used for quality control procedures in examining the accuracy and precision of the quantitative in vitro determination of Cystatin C with CysC (GFR) test kit on smart or CUBE laboratory photometer. The CysC (GFR) controls are buffer based and are supplied in a liquid, ready-to-use form.

Storage Instructions and Stability of the Control Kit

Unopened vials are stable at 2 - 8° C until the expiration date on the label. Thoroughly mix the content of the vial before each use by gently inverting it. Once opened vials are stable until the expiry date when stored tightly capped at 2 - 8° C.

Dropper bottles:

We recommend not to pipette directly from the original bottle to avoid contamination.

Use 3 - 5 drops into a Eppendorf or similar cup and pipette out of this cup.

Dispose the cup after use.

Round bottles:

Use only disposable pipette Tips to avoid contamination.

Preparation

Mix carefully before each use.

Controls are liquid and ready-to-use.

Waste Management

Please refer to your local legal requirements.

Warnings and Precautions

This control kit is for in vitro diagnostic use only.

Product contains highly purified bovine source material from non-BSE countries. Possible infectious agents in the material have been inactivated. Because no method can offer complete assurance as to the absence of infectious agents, this material should be handled as though capable of transmitting infectious disease.

Do NOT INGEST.
Avoid contact with skin and eyes.
Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds.

Take the necessary precautions for the use of laboratory reagents.

Compare the results with the values indicated in the lot specific value sheet included in the control kit.

If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor. Do not perform patient samples if the controls are outside accepted limits.

Check performance of the test according to national requirements of external interlaboratory tests (EQUAS).

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


Processing of CysC (GFR) control measurement

1.

1.1 **smart**  **1.2**  **1.3**  **1.4** **smart**  **CUBE smart** 


2.


2.1 **smart** **Sex:**
 male child
 female Control

CUBE smart Name: 6 ID: 6
Sex: Control
Samplertype: Male

2.2 **smart** **Probenart:**
 Serum
 Blut

CUBE smart Name: 29 ID: 29
Sex: Control
Samplertype: Serum

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English

ATTENTION!

Allow single test min. 10 minutes to warm up to room temperature before use!

1. Preparation of the test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette into test kit rack
- 1.3 Place ERS cap into test kit rack
- 1.4 Press measurement button, enter required information using the touch screen

2. Control preparation

- 2.1 Sex-Menu: Select "Control"
- 2.2 Sampletype-Menu: Select "Serum"
- 2.2 Dispense some drops of control liquid into cup and aspirate 20 µl using the pipette
- 2.3 Dispense 20 µl control liquid into ERS cuvette INTO LIQUID.

3. Control processing

- 3.1 Apply ERS cap firmly onto ERS cuvette
- 3.2 Place ERS Cartridge into laboratory photometer
- 3.3 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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Technical details subject to change without notice.

Doc-ID: bdce043d-94c6-45e8-ac7e-58bad6a86a0d ; latest revision: DP1.0 ; date of print: 2013-08-07