HemoCue® 201 DM Docking Station Manual

EN	HemoCue® 201 DM Docking Station Manual	4
CS	HemoCue® 201 DM Příručka k dokovací stanici	12
DE	HemoCue® 201 DM Docking Station – Bedienungsanleitung	20
DA	HemoCue® 201 DM Manual til dockingstation	28
ES	HemoCue® 201 DM Manual de la estación de conexión	36
FI	HemoCue® 201 DM Telakointiaseman käyttöohjeet	44
FR	HemoCue® 201 DM Manuel de la station de transfert	52
LT	HemoCue® 201 DM Prijungimo stotelės vadovas	60
NL	Handleiding voor HemoCue® 201 DM Docking Station	68
NO	Bruksanvisning for HemoCue® 201 DM Docking Station	76
PL	HemoCue® 201 DM Instrukcja obsługi stacji dokującej	84
SV	Bruksanvisning till HemoCue® 201 DM dockningsstation	92
SL	Priločnik za uporabo priključne postaje HemoCue® 201 DM	100
TR	HemoCue® 201 DM Yerleşirme İstasyonu Kılavuzu	108
EL	Εγχειρίδιο του HemoCue® 201 DM Docking Station	116

Docking Station Overview

Front

The Docking Station consists of the following items:

- Tracks (1) for the Analyzer
- Power (2) for the Analyzer
- LED (3)


Materials provided

- HemoCue 201 DM Docking Station (Primary or Secondary)
- HemoCue 201 DM Docking Station Manual (this document)

Back/Side panel – Primary Docking Station

The following items are located on the back/side panel:


- Power inlet (1) for the AC Adapter
- USB port (2) for setting up the Primary Docking Station
(Different from the USB port in the Secondary Docking Station.)
- Network port (3) for data communication via a LAN.

Back/Side panel – Secondary Docking Station

The following items are located on the back/side panel:


- Power inlet (1) for the AC Adapter
- USB port (2) for data communication directly with the PC
(Different from the USB port in the Primary Docking Station.)

Placing the Analyzer in the Docking Station

Always slide the Analyzer into and out of the Docking Station by means of the Tracks (1). Make sure the Analyzer is fully inserted.

Never try to lift the Analyzer out of the Docking Station or press the Analyzer downwards into the Docking Station. This may damage the casing and power outlets of the Analyzer and/or the Docking Station.


Docking Station General

For information regarding use of a HemoCue 201 DM Analyzer, refer to Instructions for Use for HemoCue Hb 201 DM Analyzer, HemoCue Glucose 201 DM Analyzer or HemoCue Glucose 201 DM RT Analyzer.

For information regarding installation of a HemoCue 201 DM System, refer to Reference Manual for HemoCue 201 DM Analyzer.


A Docking Station not connected to other Docking Stations, is referred to as a *Single Docking Station*.

Up to five Docking Stations can be connected to each other. These are referred to as *Multiple Docking Stations*.

Only the Docking Station positioned furthest to the left can be connected to an AC Adapter. Power is supplied to the other docked analyzers through this connection.

Note! Primary and Secondary Docking Stations are two physically different products.

- 1 Primary or Secondary Docking Station
- 2 Secondary Docking Station

A Primary Docking Station can be connected to a LAN and both receive/send data itself and pass data from/to a Secondary Docking Station that is connected to the Primary. When a Primary Docking Station is connected via a USB connection, only the TCP/IP settings can be changed. To change settings use PDS Configuration tool or HemoCue DMS. No other information can be sent to/from the Analyzer.

A single Secondary Docking Station cannot be connected to a LAN. In order to communicate via a LAN it must be connected to a Primary Docking Station. A Secondary Docking Station can communicate directly to a PC via a USB connection.

Both Primary and Secondary Docking Stations can be used stand alone.

A set of Multiple Docking Stations, (maximum five), consists of either one Primary and a number of Secondary Docking Stations or only Secondary Docking Stations. If a Primary is included, it must be positioned furthest to the left (1) (when facing the front of the Docking Station).

Two Primary Docking Stations cannot be connected together.

 The analyzer (including the docking station) shall only be connected to heavily restricted networks. Do not connect the analyzer or docking station to an untrusted network, application or device. Do not leave the analyzer and docking station exposed or unattended.

Getting started – Single Docking Station


Connect to a PC or a LAN

Note! This procedure should be performed by a network administrator. Important! Always keep electrical equipment turned off during connection procedures.

A Single Primary Docking Station can be connected to a LAN or to a PC via the USB.

Note! When a Primary Docking Station is connected via a USB connection, only the TCP/IP settings can be changed. To change settings use PDS Configuration tool or HemoCue DMS. No information can be sent to/from an analyzer docked in the Primary Docking Station when the Primary Docking Station is connected to a PC via USB.


A Single Secondary Docking Station can be connected directly to a PC via the USB, but not to a LAN.

Connect Primary Docking Station to PC

To change settings on a Primary Docking Station, see section Docking Station Settings in the Reference Manual for the HemoCue 201 DM Systems. The Primary Docking Station must be connected to a PC directly via the USB port.

- Insert the USB cable into the USB port (1).
 - Insert the USB cable into the PC's USB port
 - Connect the AC Adapter as described in section *Power Source*.
- Note! The AC adapter must be connected after task a) and b).


Connect a Primary Docking Station to a LAN

A Primary Docking Station must be connected to a LAN for data communication. The PC hosting the DMS Software, or an Observation Reviewer must also be connected.

Docking Stations labelled PDS+ and SDS+ supports network communication speed 10/100Mbps.


- Insert one end of the Network cable into the Network port on the Docking Station (2).
 - Insert the other end of the Network cable into a Network port on the LAN (not directly into the PC).
 - Connect the AC Adapter as described in section *Power Source*.
- Note! The AC adapter must be connected after task a) and b).

Connect a Secondary Docking Station to a PC

A Single Secondary Docking Station must be connected via the USB to the PC hosting the DMS Software for data communication. No settings can be changed in the Secondary Docking Station.

- Insert one end of the USB cable into the USB port of the Docking Station (1).
- Insert the other end of the USB cable into the PC's USB port.

Power source

Connecting the AC Adapter

! Only use the AC Adapter that is recommended by HemoCue.

- a) Insert the AC Adapter's DC plug (2) into the Power inlet (3).
- b) Plug the AC Adapter (1) into an electrical outlet.

The Docking Station (Primary or Secondary) furthest to the left is connected to the Power source.

The rest of the Docking Stations in a multiple set are connected to the Power source via this Docking Station and cannot be connected separately.

Getting started – Multiple Docking Stations

Connecting Docking Stations

Important! Docking Stations should be connected to each other before connecting to the Power source, the PC or the LAN.

Multiple Docking Stations, (max 5), consists of either one Primary and a number of Secondary Docking Stations or a number of Secondary Docking Stations only.

Note! If a Primary Docking Station is included, it must be positioned furthest to the left (when facing Docking Station fronts).

Primary Docking Stations labelled PDS⁺ shall only be connected to a Secondary Docking Station labelled SDS⁺. Secondary Docking Stations labelled SDS⁺ shall only be connected to Docking Stations labelled PDS⁺ or SDS⁺. Note! Malfunction will occur if non-labelled Docking Stations are connected to PDS⁺ or SDS⁺.

Note! Two Primary Docking Stations cannot be connected together.

Note! Max 4 Secondary Docking Stations should be connected to a Primary Docking Station.


- a) Remove Side Trim (2) from Primary Docking Station by removing the screw* from the bottom side (if present), and then gently pressing Side Trim inwards, while lifting it upwards. Male connector (1) is uncovered.

- b) Remove the Bottom Trim (4) from the Secondary Docking Station to uncover the female connector (3). This is done by gently inserting a small screwdriver** into the small openings and gently lifting the trim upwards.

*Use a Torx screwdriver, size T10.


**Use a flathead screwdriver.


c) Gently press the Secondary Docking Station (5) downwards onto the Primary Docking Station (6), fitting the Female Connector and the Male connector together.

Note! The Docking Stations must be securely connected to each other in order to ensure good electrical contact.

d) If more Secondary Docking Stations are to be added, remove the Side Trim from the previously added Docking Station and the bottom trim of the new Secondary Docking Station, follow step a) and b). Repeat the procedure in step c) with the new Secondary Docking Station.

Important! If not connected to another Docking Station the Side Trim (2) on the Docking Station should always be used to protect the Connector (1).

Docking Station LED

Steady green light

A steady green light from the LED (1) indicates that the Docking Station (3) is receiving power and that the Analyzer Battery is fully charged.

Flashing green light

A flashing green light from the LED (1) indicates that the Analyzer Battery in the docked Analyzer (2) is charging.

Steady red light

A steady red light from the LED (4) indicates an internal error within the Docking Station. See section Troubleshooting.

Note! Not valid for unlabelled Secondary Docking Stations (non-SDS⁺).

Flashing red light

A flashing red light from the LED (4) indicates an external communication error. See section Troubleshooting.

Note! Not valid for unlabelled Secondary Docking Stations, (non-SDS⁺).

Troubleshooting Guide

Symptom	Explanation	Action
No transfer of data	<ol style="list-style-type: none"> 1) No network communication with OR 2) No network communication with PC 3) Communication error 4) No communication 5) No communication between the first and second Docking Station 6) No communication between the Secondary Docking Stations 7) No communication with any Docking Station 	<p>If multiple Docking Stations are used, make sure not to mix Docking Stations labelled PDS⁺/SDS⁺ with unlabelled Docking Stations.</p> <ol style="list-style-type: none"> 1a) Check the network parameters. 1b) Check that the LED on the Docking Station does not show a red light (stable or flashing). 2a) Check the network parameters. 2b) Check that the LED on the Docking Station does not show a red light (stable or flashing). 3) More than 5 Docking Stations are docked together. 4) Check that a USB connection is not used for communication in a Primary Docking Station. 5) Check that two Primary Stations are not docked together. 6) Check the electrical connections. The LED on the docking station should give either a stable green light or a flashing green light. 7) Check the electrical connection. The LED on the docking station should give either a stable green light or a flashing green light. <p>If problems remain contact distributor.</p>
No transfer of data via USB	No USB-communication	<ol style="list-style-type: none"> a) Check that the Analyzer is properly docked. b) Check the cable connections in both the Docking Station and the PC. c) Check that the Docking Station is not a Primary Docking station. d) Check that all Docking Stations that should communicate are connected with a USB cable to the PC. e) Check that the DMS Software is properly installed and that the USB driver is installed. f) Check the electrical connection. The LED on the docking station should give a stable or flashing green light.
Analyzer not charged	No charging of the Battery	<ol style="list-style-type: none"> a) Check that the Analyzer is properly docked. b) Check that the green LED on the Docking Station gives a flashing green light when docking the Analyzer. c) Replace the Battery.

Technical Specifications

Dimensions	206 × 135 × 61 mm (8.10 × 5.30 × 2.40 inches)
Weight	566 g (1.24 pounds)
System operating temperature	See relevant Instructions for Use
Docking station storage and transport humidity and temperature	< 90% RH (non-condensing) 0 - 50 °C (32 - 122 °F)
Docking station operating humidity	< 90% RH (non-condensing) 15 - 25 °C (59 - 78 °F) < 75% RH (non-condensing) 26 - 30 °C (79 - 86 °F)
Electrical Rating	Primary Docking Station: 6 V, 2.5 A Secondary Docking Station: 6 V, 0.5 A
AC adapter Only use adapters recommended by HemoCue.	FW7556M/06: Input: 100V~ – 240 V AC/50–60 Hz/400-200 mA Output: 6 V DC/2.5 A FW8001/06: Input: 100V~ – 240 V AC/50–60 Hz/400-200 mA Output: 5.9 V DC/3.0 A
Pollution degree	2
Over voltage category	II
Atmospheric pressure	700 hPa to 1060 hPa
Equipment not suitable for use in the presence of flammable mixtures The HemoCue 201 DM Docking Station is tested according to IEC 61010-1, IEC 61326-1 and IEC 61326-2-6	

Electromagnetic Compatibility

The HemoCue 201 DM Docking Stations complies with emission and basic immunity requirements described in IEC 61326-1.

The HemoCue 201 DM system complies with the emission and immunity requirements described in IEC 61326-2-6. When using the analyzer standalone or together with one Primary Docking Station it has been designed and tested to CISPR 11 Class B. When using the analyzer with a standalone Secondary Docking Station or one or several Secondary Docking Stations connected to a Primary Docking Station it has been designed and tested to CISPR 11 Class A. In a domestic environment it may cause radio interference, in which case, you may need to take measures to mitigate the interference.

The electromagnetic environment in which the HemoCue 201 DM systems will be operated should be evaluated prior to operation of the device. Do not use the HemoCue 201 DM systems in close proximity to sources of strong electromagnetic radiation (e.g. mobile phones and other wireless communication devices) as these can interfere with the proper operation.

Note: It is the manufacturer's responsibility to provide equipment electromagnetic compatibility information to the customer or user.

Note: It is the user's responsibility to ensure that a compatible electromagnetic environment for the equipment can be maintained in order that the device will perform as intended.

Warranty

The HemoCue 201 DM Docking Station carries a 24-month warranty from the day of receipt. After the warranty period, service/repair is carried out at fixed prices. Any other use of the System than recommended by the manufacturer, will void the warranty.

Service and Disposal

The 201 DM Docking Station should be cleaned as recommended under Maintenance in the Instructions for Use for HemoCue Hb 201 DM Analyzer, HemoCue Glucose 201 DM Analyzer or HemoCue Glucose 201 DM RT Analyzer, respectively, prior to service or disposal. Consult your local environmental authorities for proper disposal.

Symbols used

Caution

Consult Instructions For Use

In Vitro diagnostic medical device


Manufacturer

Direct Current

Efficiency level

Ethernet

USB

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

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HemoCue AB
HemoCue AB | Kuvettgatan 1 | SE-262 71 Ängelholm | Sweden