## Instructions for Use – PA-AST Panel U-0501

### Information

Please consult the full version of the *Instruction for Use* (CG396173) of the PA-100 AST System before use.

## Intended use of the PA-AST Panel U-0501

PA-AST Panel U-0501 together with the PA-100 Analyzer is intended for automated *in vitro* qualitative presumptive determination of bacteriuria and qualitative Antibiotic Susceptibility Testing (AST) of most common bacterial pathogens in Urinary Tract Infections (UTI), directly from urine sample. The results are intended for use as an aid in the diagnosis and selection of antibiotic treatment of UTI in conjunction with other clinical signs and symptoms. The system comprising the PA-AST Panel U-0501 together with the PA-100 Analyzer is for *in vitro* diagnostic use by Health Care professionals in Near Patient Testing. The PA-AST Panel U-0501 is only to be used with the PA-100 Analyzer.

## Indications for use

Bacteriuria detection and AST performance of the system has been evaluated with the most common uropathogenic bacterial species listed below:

**Gram negatives:** *Escherichia coli*, *Klebsiella pneumoniae*, *Proteus mirabilis*


**Gram positives:** *Staphylococcus saprophyticus*, *Enterococcus faecalis*


The PA-AST Panel U-0501 is indicated for *in vitro* presumptive detection of bacteriuria (reported as **positive bacteriuria**), that is the presence of actively growing bacteria in patient urine samples at a concentration above the cut-off value of  $5 \times 10^4$  Colony-Forming Units/mL (CFU/mL). The absence of actively growing bacteria will result in **negative bacteriuria**. Note that reported negative bacteriuria does not mean that there are no bacteria in the urine. For example, concentrations of actively growing bacteria below the cut-off can result in negative bacteriuria.

The PA-AST Panel U-0501 is indicated for *in vitro* Antimicrobial Susceptibility Testing (AST) for the following five antibiotics for urine samples that are reported positive for bacteriuria by the PA-100 AST System:

- Amoxicillin-Clavulanic acid (Amoxi-Clav)
- Ciprofloxacin
- Fosfomycin
- Nitrofurantoin
- Trimethoprim

## Symbols

Symbol	Title of Symbol
	CE mark
	In vitro diagnostic medical device
	Catalogue number
	Batch code
	Manufacturer
	Do not use if package is damaged and consult <i>instructions for use</i>

Symbol	Title of Symbol
	Temperature limit
	Consult instructions for use
	Use-by date
	Do not re-use
	Contains sufficient for 25 tests

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