QUICK REFERENCE GUIDE
Easy operation – fast results
ENGLISH - DEUTSCH - FRANÇAIS - ITALIANO - ESPAÑOL

Index

English
4-11

Deutsch
12-19

Français
20-27

Italiano
28-35

Español
36-43

English

Deutsch


Français

Italiano

Español

Unpacking and Setup

- 1 Remove the *Piccolo xpress* chemistry analyzer from its shipping carton.
- 2 Place the analyzer on a level surface that is free of hair, dust, and other contaminants. Do not place the analyzer near a sunny window or any other heat source. Make sure the analyzer is at least 20 cm away from any wall to provide adequate ventilation and access to the power connection, Ethernet and USB ports.
- 3 Check the components received with the *Piccolo xpress* against the following figure to make sure you have everything required to set up the analyzer.

Note: The shipping carton also includes following items:
USB cable, Abaxis Driver CD, Multilingual Operator's Manual CD, warranty card, and accessory kit.

- 4 Complete the warranty card and mail it to Abaxis within 10 days to start the warranty period.
- 5 Open the printer cover.
- 6 Unwrap the paper roll, then unroll several centimeters.
- 7 Put the paper into the printer so that it unrolls from the bottom towards the front of the analyzer as shown. Make sure several centimeters of paper extend out of the printer slot.
- 8 Press the cover closed until it locks into place.
- 9 Pull gently on the end of the paper until it is taut.

- 10 Plug the power jack into the analyzer, then plug the detachable power cord into the power adapter and into a grounded electrical outlet.

Note: *To prevent power surges or drain, do not plug the analyzer into the same circuit as a centrifuge or any other high-current device. Abaxis also recommends using a suitable surge protector.*

- 11 Press the power button to turn on the analyzer. The display then shows the following:


- 12 After passing the self-test and reaching operating temperature, the analyzer is ready to run the first reagent disc, and displays **Analyze**. See *“Quality Control Features”* on page 6-1 of the operator’s manual for details about the analyzer self-test. The analyzer may require additional time to reach operating temperature.
- 13 Check the analyzer date and time to ensure they are correct. See *“Changing Date and Time”* on page 4-15 of the operator’s manual for directions.
- 14 Reference ranges are already preset by the factory. The range values can be changed using the Customizing Reference Range feature described in *“Customizing Reference Ranges”* on page 4-2 of the operator’s manual.

Please make sure to be familiar with the system before running any patient samples.

Filling the Reagent Disc

Using the provided 100 μL micropipette.


- 1 Take a reagent disc out of its pouch. Discs can be used directly from the refrigerator without warming.
- 2 Attach a new tip to the end of the pipette.
- 3 Fill the pipette without air bubbles or air gaps with sample.
- 4 Tilt the disc to 45° with the sample port above the fill line, so that the entire sample will flow into the sample chamber.
- 5 Place the pipette tip into the disc's sample chamber and release the sample. Make sure to keep the pipette plunger pressed down until the pipette tip is removed from the sample port.
- 6 Discard the pipette tip into a biohazard container.
- 7 Clean the reagent disc. Use a lint-free tissue to remove any sample spilled on the outside of the disc, taking care that the tissue does not withdraw any sample from the sample port. Dispose tissue in a biohazard container.

Beside the micropipette, any other transfer device can be used to dispense approximately 100 μL of sample into the disc via the sample chamber.

Note: *When using lithium heparin whole blood tubes, gently invert them several times before use.*

Running a Patient Sample

- 1 Turn on the analyzer by pressing the Power button on the front of the analyzer. The analyzer starts up, then performs a self-test. During the warm up time the display shows **Warning**. When the analyzer reaches operating temperature, it displays **Analyze** on the Home screen.
- 2 Select **Analyze** on the touchscreen to open the disc drawer.
- 3 Place the disc in the recessed area in the drawer.

- 4 Select **Close** on the touchscreen to close the disc drawer.
- 5 Select **Patient** or **Control** to run a patient sample or serum control.
- 6 Enter an ID number for the sample (up to 14 characters), then select **Done**.
- 7 The analyzer then checks the disc type, and begins processing the sample.
- 8 When the sample is finished processing, the analyzer stores the results and shows that the analysis is complete.
- 9 By default, the analyzer automatically prints the results of the analysis.
- 10 Select **Open** to open the disc drawer and discard the used reagent disc into a biohazard container.
- 11 To analyze another sample, insert a new reagent disc and repeat the above procedure or when finished, select **Close** to close the disc drawer and to get the analyzer returned to the home screen.

Note: *The Piccolo xpress chemistry analyzer accepts lithium-heparinized whole blood, plasma, or serum for patient samples. Lithium heparin is the only anticoagulant recommended for use with the Piccolo xpress.*

Interpreting Results

English

piccolo xpress			
MetLyte Plus CRP			
02 Aug 2013	02:44 AM		
Sample Type:	Patient		
Patient ID:	020813		
Disc Lot Number:	3192BA2		
Serial Number:	0000P20026		
.....			
GLU	111	73-118	mg/dL
BUN	16	7-22	mg/dL
CRE	0.5 *	0.6-1.2	mg/dL
CK	79	30-380	U/L
NA+	139	128-145	mmo1/L
K+	4.5	3.6-5.1	mmo1/L
CL-	100	98-108	mmo1/L
tCO2	26	18-33	mmo1/L
CRP	< 5.0	0.0-7.5	mg/L
.....			
QC	OK		
HEM	0	LIP	0
		ICT	0

Head of result printout:

- Disc type
- Date/ Time
- Sample Type
- Sample-ID
- Lot#
- Serial#

Test results:

- Analyte
- Analyte Concentration
- Reference range
- Units

QC OK:

iQC™ checks passed

HEM/LIP/ICT:

Qualitative indices for Sample interferences

Symbols printed on the printout

- * result outside the reference range
- < result below the dynamic range of the analyzer
- > result above the dynamic range of the analyzer
- c calculated value

QC OK all iQC™ checks passed successfully

Qualitative indices indicate the degree of hemolysis (HEM), icterus (ICT), and lipemia (LIP) found in the sample on a scale of: 0 (clear), 1+ (slight), 2+ (moderate), and 3+ (gross).

The analyzer features an internal quality control (iQC™) system that suppresses values and displays warning or error messages when a problem occurs. The iQC™ conducts several quality and sufficiency checks with every run while monitoring all technical components of the analyzer. A suppression or warning does not necessarily mean that your instrument has a technical fault. It proves that your iQC™ is working properly and ensures that only accurate results will be released.

- The symbols ~~~ are printed in place of numbers when a result cannot be determined. The result may be suppressed due to improper mixing of a reagent bead with diluted sample, a nonlinear reaction, an endpoint of a particular reaction not reached, a concentration outside the analyzer's capabilities or thermal stress. When a chemistry is suppressed (~~~), the analyzer prints an error report.
- HEM, LIP, or ICT is printed in place of the analyte concentration if hemolysis, lipemia or icterus have adversely affected the results. LIP is also printed if both lipemia and icterus have affected a particular analyte. HEM is also printed if hemolysis and icterus, hemolysis and lipemia, or hemolysis, lipemia, and icterus have affected a particular analyte. Examine the sample indices to determine if more than one interferent is affecting a particular result.
- CONFIRM LOW RECOVERIES with an exclamation point (!) next to every analyte indicates that at least one analyte has a lower concentration than normally expected. If this occurs, re-run the sample – if the message persists, the sample may be problematic.
- Warnings and error messages include an internal error code that will assist the Abaxis Technical Support in diagnosing the issue. Record the error message and/or print the error report (see chapter 5.2 of the operator's manual) before contacting Abaxis Technical Support:

techsupport@abaxis.de

Cleaning the Analyzer

Spills

Observe universal precautions when cleaning spills on the analyzer. Use a 10% bleach solution (1 part bleach plus 9 parts water) to clean spills, following the standard cleaning guidelines in “Cleaning the Analyzer” on page 9-1 of the operator’s manual.

Display

Clean the analyzer’s screen periodically using a soft, lint-free cloth dampened with a glass-cleaning fluid or window cleaner.

The screen can be disinfected using a 10% bleach solution: apply the solution to a lint-free cloth, then wipe the screen.


CAUTION: Do not use any cleaner containing alcohol.
Do not spray cleaner directly onto the display – dampen the cloth instead

Air Filter

The air filter in the back of the analyzer should be cleaned at least twice per year. Check the air filter more frequently if the analyzer is located in an environment with excessive dust or dirt.

To clean the air filter:

- 1 Unplug the analyzer and remove the power cord from the back of the analyzer.
- 2 To reach the filter, remove the air filter cover at the back of the analyzer by pulling.
- 3 Wash the filter in warm soapy water and dry it completely before reassembling.
- 4 Place the clean, dry air filter into its cover. Afterwards push the 4 clips back to their position on the analyzer.
- 5 Plug the power cord back into the analyzer.

Piccolo xpress™
Panels

Specification

	General Chemistry 6	General Chemistry 13	Biochemistry Panel Plus	Basic Metabolic Panel	Basic Metabolic Panel Plus	Comprehensive Metabolic Panel	Electrolyte Panel	Metlyte 8 Panel	Metlyte Plus CRP	MetLac 12 Panel	Lipid Panel	Lipid Panel Plus	Liver Panel Plus	Hepatic Function Panel	Renal Function Panel	Kidney Check	
ALB		●	●			●				●		●	●	●			Albumin
ALP		●	●			●						●	●				Alkaline Phosphatase
ALT	●	●	●			●					●	●	●				Alanine-Aminotransferase
AMY		●	●									●					Amylase
AST	●	●	●			●					●	●	●				Aspartate-Aminotransferase
DBIL													●				Direct Bilirubin
BUN	●	●	●	●	●	●	●	●	●	●					●	●	Blood Urea Nitrogen
Ca		●	●	●	●	●				●					●		Total Calcium
CK							●	●									Creatine Kinase
Cl ⁻				●	●	●	●	●	●	●					●		Chloride
CRE	●	●	●	●	●	●	●	●	●						●	●	Creatinine
CRP			●					●									C-reactive Protein
GGT	●	●	●										●				Gamma-Glutamyl-Transferase
GLU	●	●	●	●	●	●	●	●	●	●	●				●		Glucose
K ⁺				●	●	●	●	●	●	●					●		Potassium
LAC										●							Lactate
LDH					●												Lactate-Dehydrogenase
Mg					●					●							Magnesium
Na ⁺				●	●	●	●	●	●	●					●		Sodium
PHOS										●					●		Phosphorus
TBIL		●				●							●	●			Total Bilirubin
tCO ₂				●	●	●	●	●	●	●					●		Total Bicarbonate
TP		●	●			●							●	●			Total Protein
UA		●	●														Uric Acid
CHOL									●	●							Cholesterol
CHOL/ HDL*									●	●							Cholesterol/HDL*
HDL									●	●							High Density Lipo Protein
nHDL*									●	●							Non-HDL*
LDL*									●	●							Low Density Lipo Protein*
TRIG									●	●							Triglycerides
VLDL*									●	●							Very Low Density Lipo Protein*

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ABAXIS Europe GmbH
ABAXIS Inc.