Control kit for quality control for the Ferritin parameter on the smart analyser.

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Kit size

32 Tests/Package
2 x 1 ml / decision level

Test kit preparation: Allow single test min. 15 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Ferritin control kit is intended for use in the clinical laboratory as a control for quantitative determination of Ferritin in the normal range and high range. It is prepared from a human serum and was found to be non reactive when tested with approved methods for HbsAg, anti HIV1/2 and anti HCV. However this material should be handled as though capable of transmitting infectious disease.

Preparation and Handling

The control is liquid and ready to use.
Thoroughly mix the content of the vial before each use by gently inverting it.
Turbidity shows bacterial contamination, the vial has to be discarded.
Always use single-use pipette tips!

Stability and Storage

Unopened vials are stable at 2-8°C until the expiration date on the label. Once opened vials are stable for at least 8 weeks when stored tightly capped at 2-8°C.

Waste Management

Please refer to local legal requirements.

Warnings and additional Informations

This testkit is for in vitro diagnostic use only. Do NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbin to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

- Y Set the Sex to "Control" before running the QC serum.
- Y Compare the results with the values indicated in the lot specific value sheet included in the control kit
- Y If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor.

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
Processing of Ferritin control measurement

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




English

1. Preparation of test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette in smart test kit rack
- 1.3 Place ERS cap in smart test kit rack
- 1.4 Press symbol on touch screen, required information is entered by use of the touch display of the analyser

2. Settings on the smart Analyser

- 2.1 Gender/Sex-Menu: Select "Control"
- 2.2 Instrument ready for start of analysis

3. Processing Ferritin control measurement

- 3.1 Dispense some drops control liquid into serum cup. Aspirate 50µl from serum cup.
- 3.2 Dispense 50µl sample INTO LIQUID in ERS cuvette
- 3.3 Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into smart analyser.
- 3.5 Start automatic sample processing by pressing the  start button.

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ATTENTION!

Allow single test min. 15 minutes to warm to room temperature!

Wertetabelle für smart Ferritin Control Kit Value sheet for smart Ferritin control kit

g

Control Value

Target	value (min)	value (max)

H

M

Eurolyser Diagnostica GmbH
www.eurolyser.com