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SUPER GL *speedy*

Elektrochemischer Analysator

REF 0330 000 000 00

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Gültig ab 1. Juni 2022

Hersteller / Copyright:

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Dr. Müller Gerätebau GmbH
Burgker Str. 133
01705 Freital
Germany

Service – Hotline:

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Dr. Müller Gerätebau GmbH
Burgker Str. 133
01705 Freital
Germany

Service hotline:

If you have any questions concerning our products
please contact us:

+49 (0)351 65 20 86 21

Fax: +49 (0)351 65 20 86 30

Email: support@diabetes-diagnostics.com

You will find the current version of this manual on our website

www.diabetes-diagnostics.com

A) Table of contents:

	Page
A..... Table of contents.....	2
B..... Table of figures.....	4
C..... Table of symbols.....	5
1..... Preface.....	6
1.1..... Introduction.....	6
1.2..... The SUPER GL speedy.....	6
1.2.1..... Basics.....	6
1.2.2..... Device and accessories.....	7
1.2.3..... Intended use.....	7
1.3..... Indication / Contraindication.....	8
1.4..... Manufacturers liability.....	8
2..... Safety.....	9
2.1..... Introduction.....	9
2.2..... Responsibility / Training of the operator.....	9
2.3..... General safety instructions.....	9
2.4..... Product-specific safety instructions.....	10
3..... Description of analyser.....	11
3.1..... Introduction.....	11
3.2..... Measuring principle.....	11
3.3..... <i>Calibration mathematics</i>	12
3.4..... <i>Consumables</i>	12
3.4.1..... <i>System Solution SUPER GL</i>	12
4..... Operation - Part 1.....	13
4.1..... Introduction.....	13
4.2..... Safety instructions.....	13
4.3..... Installing the device.....	14
4.3.1..... <i>Installing accessory printer TPU-80</i>	15
4.4..... Initial operation <i>of the device</i>	16
4.5..... Preparing the measuring process.....	19
4.5.1..... Basics.....	19
4.5.2..... Sample preparation.....	19
4.6..... Measuring.....	21
4.6.1..... Rotor operation.....	21
4.6.2..... Chain operation.....	21
4.6.3..... Calibration.....	22
4.6.4..... Controls.....	23
4.6.5..... Method.....	23
4.7..... Switching off the device.....	23

Table of contents

	Page
5	24
5.1	24
5.2	24
5.3	24
5.3.1	24
5.3.2	24
5.3.2.1	24
5.3.2.2	25
5.3.2.3	26
5.3.2.4	26
5.3.2.5	26
5.3.3	26
5.3.3.1	26
5.3.3.2	27
5.3.3.3	28
5.3.3.4	28
5.3.3.5	28
5.3.3.6	28
5.3.3	29
5.3.4	29
6	30
6.1	30
6.2	30
6.3	30
6.3.1	30
6.3.2	31
6.3.3	32
6.3.4	33
6.3.5	33
6.4	34
6.4.1	34
6.4.2	34
6.4.3	35
6.4.4	36
7	37

B) Table of figures

		Page
Overall view SUPER GL speedy.....	Fig. 1.1.....	6
<i>Package content</i>	<i>Fig. 1.2</i>	7
<i>Optional accessories</i>	<i>Fig. 1.3</i>	7
Flow schema.....	Fig. 3.1.....	11
<i>Reactions in the sensor for glucose</i>	<i>Fig. 3.2</i>	11
<i>Reactions in the sensor for lactate</i>	<i>Fig. 3.3</i>	11
Connections SUPER GL speedy.....	Fig. 4.1.....	14
<i>TPU-80 printer</i>	<i>Fig. 4.2</i>	15
<i>Back side view TPU-80</i>	<i>Fig. 4.3</i>	15
Touch Screen.....	Fig. 4.4.....	16
Main program branches.....	Fig. 4.5.....	18
Sample preparation using open-end capillary.....	Fig. 4.6.....	20
Sensor holder closed.....	Fig. 6.1.....	31
Sensor holder open.....	Fig. 6.2.....	31
View supply and waste container.....	Fig. 6.3.....	32
<i>Table of displayed and printed warnings</i>	<i>Fig. 6.4</i>	34
Table of errors during measuring.....	Fig. 6.5.....	34
List of technical data.....	Fig. 7.1.....	37

C) Table of symbols

General used symbols

Symbol	Description
	<i>Biohazard</i>
	Follow instruction manual
IVD	In-vitro-Diagnostica
	Manufacturer
	<i>Distributor</i>
	<i>Date of manufacture</i>
	<i>For near-patient testing</i>
CE	CE Compliance
SN	Serial number
REF	Item number
Cont.	Contents of package
LOT	Batch number
	Use before
	Storage temperature

Symbols in instruction manual

Symbol	Description
	Attention
Bold	Very important notes
italics	<i>Modifications</i>

1 Preface

1.1 Introduction

Congratulations on purchasing the SUPER GL speedy analyser. We hope you will find working with your analyser satisfying and successful.

In the following chapter "The SUPER GL speedy" you will find a first overview of your analyser: what parameters you can measure, what further devices and accessories belong to your analyser, and an overview of the device's functionality.

Furthermore, you will receive information on safety, on liability and warranty, on indications *and* contraindications of your analyser.

For further and more detailed information, please read the corresponding chapters.

1.2 The SUPER GL speedy

1.2.1 Basics

The SUPER GL speedy was designed using the latest technology along with decades of experience in the area of production of clinical-chemical analysers.

Fig. 1.1 Overall view SUPER GL speedy

It fulfils all legal specifications with regard to design and production that are required of all devices used in clinical laboratories. The compliance with the valid norms and statutes is documented through the visibly attached CE-Label. The CE-Label signifies compliance with all pertaining laws and regulations and consequently safety and confidence.

Please refer to the corresponding Declaration of Conformity.

By employing an altogether newly developed technology for the determination of glucose and lactate, it is possible to fulfil all requirements of quality assurance (e.g. RiliBäk (Guidelines of the Federal General Medical Council for Quality Assurance in Medical Laboratories)) in medical laboratories while maintaining easy handling and minimum operating effort. All users are thus able to achieve analysis results that meet the quality demands.

1.2.2 Device and accessories

Deliverables:

Designation	Quantity
SUPER GL speedy	1
Sample rotor	1
Waste container	1
Power connection cable	1
Power supply adapter for analyser	1
Accessory printer TPU-80 including printer cable	1
Instruction manual	1

Fig. 1.2 Package content

Optional

EDP cable (SUB-D9)	1
EDP cable (USB)	1
Sample rotor 2	1

Fig. 1.3 Optional accessories

1.2.3 Intended use

The SUPER GL speedy is a fully automated analyser for professional biochemical in-vitro diagnostics in clinical laboratories and laboratory-like environments. Electrochemical biosensors, that work according to the enzymatic-amperometric principle can be used. The device must only be used and operated by trained personnel. Operating the device for self-testing is expressly forbidden.

1.3 Indication / Contraindication

Indication:

The SUPER GL speedy analyser is an automatic analyser for the determination of *different biochemical parameters* in prediluted samples, e.g. in hemolysed blood samples. *Actually, biosensors for the quantitative measurement of glucose and/or lactate are available. Therefore the following chapters are focused on this two parameters.*

For the selection of sample material, please refer to package inserts of the biosensors and sample cups. Use only the specified sample materials, else faulty measured results may occur.

The sample material is taken out of closed reaction cups that are placed on a sample rotor. The containers for the system solution and waste solution are located inside the device.

The device has the following features:

- automatic *quick* sample measurement (*STAT*) or sample series measurement up to 60 samples
- automatic detecting the barcode
- automatic calibration *modes*
- serial *RS 232* printer interface
- serial *RS 232* EDP interface

Contraindication:

Using unsuitable sample material can result in faulty measuring results. If in doubt, call the manufacturer!

Operating the device for *self*-testing is expressly forbidden!

1.4 Manufacturers liability

Legal liability and warranty claims are expressly excluded in the following cases:

- gross negligence or wilful damage of the device, parts thereof or *consumables*
- unauthorised opening of the device by untrained personnel (without proper service training)
- force majeure (e.g. stroke of lightning, water damage, fire)
- nonobservance of user manual and package inserts

2 Safety

2.1 Introduction

The following chapters concern the safety of the person operating the device.

Read these chapters carefully **PRIOR** to starting up the device because they contain general safety warnings, warnings concerning the personal safety of the person operating the device, and warnings for the protection of the device.

Displaying the following safety warnings does not release the person operating the device from adhering to the safety measures of the facility.

2.2 Responsibility / Training of the operator

- The SUPER GL speedy analyser must only be used and operated by trained personnel. *Trained personnel are people with an education with medical or laboratory scope and who have gained expert knowledge by studying this instruction manual and hands on at the SUPER GL speedy. A training course by an employee of the manufacturer or of an authorised distributor can support gaining knowledge.*
- Every user is responsible for adhering to safety, health and legal regulations, and operating the device only according to its intended use.
- Interpreting the results and diagnosing on that basis must be left to a medical practitioner. Operating the device for *self*-testing is expressly forbidden.

2.3 General safety instructions

- Prior to operating the device, read the entire instruction manual especially the instructions for sample taking. If you have any questions, please contact the manufacturer or authorised distributor.
- Every person working with the device must be acquainted with the relevant safety rules prior to operating the device and these rules should be kept at hand all the time.
- Please pay attention to all general safety rules for laboratories such as wearing protective gloves, and the applicable disinfection and hygiene regulations.
- *If there is any accidental skin contact with potentially infectious substances, e.g. human samples, disinfect the affected parts of the body with a suitable disinfectant solution.*
- To avoid risk of electric shock, do not place the device or power supply in water or other liquids! If the cable or the power supply adapter is damaged in any way, you must not continue using the power supply. Never touch the plug of the power supply adapter with wet hands. The power supply adapter must only be used indoors and must be protected from humidity.

2.4 Product-specific safety instructions

- The device may only be used for the intended use with special attention to the defined usage restrictions and constraints that have to be strictly adhered to (if need be, contact the manufacturer).
- Operate the device only on smooth, horizontal surfaces. Avoid variations in temperature, drafts, direct sun light, and vibrations. These can result in faulty measuring values.
- In case of malfunctions, stop operating the device immediately! Prior to continuing to operate the device, read the notes concerning cleaning, error messages and troubleshooting. After consulting the manufacturer or authorised distributor you may ship the device for repairs to the manufacturer or authorised distributor.
- Use only original accessories and spare parts to avoid damage to device and people. Repairs must only be conducted by the manufacturer or by companies authorised by the manufacturer!
- The use of reagents and *consumables* that are not expressly recommended by the manufacturer can cause severe measuring errors and malfunctions and is therefore not permissible.
- If the user opens the device without authorisation, the user shall not be entitled to any rights concerning the liability for the device and damages caused thereby.

3 Description of analyser

3.1 Introduction

This chapter describes the analyser's measuring principle, *the calibration mathematics and the consumables*. For detailed instructions and descriptions of the device, please refer to chapter Operation.

3.2 Measuring principle

The determination of *analytes* with the SUPER GL speedy is based on an electro-chemical measuring principle with a biosensor. With the help of an analyser pump system solution, calibration solution, control material or patient's material is conveyed through a sensor. The electrodes inside the sensor are separated from the liquid by semipermeable layers containing *e.g.* immobilized enzymes. The following figure shows the flow schema. *The electrochemical reaction takes place inside the biosensor.*


Fig. 3.1 Flow schema

The next two figures show these reactions as example for the measurement of glucose and lactate.


Fig. 3.2 Reactions in the sensor for glucose


Fig. 3.3 Reactions in the sensor for lactate

3.3 Calibration mathematics

The change of the electrical current on the biosensor will be used to calculate the concentration of the analyte based on a calibration curve. Because the correlation between concentration of analyte and measuring signal is linear a linear calibration mathematics is used. The concentration of the sample will be calculated by the following equation:

$$C_{\text{Sample}} = \frac{C_{\text{Calibrator}}(I_{\text{Sample}} - I_{\text{SL}})}{I_{\text{Calibrator}} - I_{\text{SL}}}$$

C_{Sample} = Concentration of sample

$C_{\text{Calibrator}}$ = Concentration of calibrator

I_{Sample} = Current of sample

$I_{\text{Calibrator}}$ = Current of calibrator

I_{SL} = Current of System Solution

3.4 Consumables

To operate the analyser following consumables are needed:

- Glucocapil reaction cups prefilled with Hemolysate-System-Solution
- capillaries or pipette for dosing sample
- Glucocapil calibration solution
- System Solution SUPER GL in 1.2 litre containers
- biosensor for the analyte
- control material

For detailed instructions on the use of these consumables, please refer to the following chapters of this manual respectively the separate product information of these consumables.

3.4.1 System Solution SUPER GL

The System Solution SUPER GL is intended for use as washing solution for automatic filling and washing the liquid system of the analyser SUPER GL and to provide a stable working environment for the biosensors.

The System Solution is not corrosive, poisonous or seriously harmful to skin. If there is any skin contact with the System Solution, rinsing with water is sufficient.

The System Solution SUPER GL must be stored protected from light and at temperatures of +2°C and above. It is ready to use and stable until the indicated expiry if protected from contamination, proper storage and handling are observed.

Containers with System Solution SUPER GL and for waste will be installed in the device. Please refer chapter 6.3.3 for installation and exchange! Empty containers should be used as waste containers.

4 Operation - Part 1

4.1 Introduction

In this part of the instruction manual, all information is included that is useful for the day-to-day operation of the device.

In part 2, all additional information is included that is important for understanding the functions, complementary functions and certain sources of possible problems.

The qualified personnel for the device must be familiar with both parts and must also have the medical knowledge to be able to interpret the acquired values correctly. Conclusions for *diagnosis and* therapy may only be drawn by a medical practitioner.

4.2 Safety instructions

As mentioned before, certain safety warnings must be heeded when operating the device to guarantee correct and faultless operation:

- The device must only be used for the described *intended use and with observance of* indication. *It* must only be used and operated by trained personnel.
- Every user is responsible for adhering to safety, health and legal regulations, and operating the device only according to its intended use.
- Interpreting the results and diagnosing on that basis must be left to a medical practitioner. Operating the device for *self-testing* is expressly forbidden.
- In daily operation, regular checks of the results should be made; if needed, an additional control measurement should be carried out.
- Do not switch off the device or disconnect it from the power supply while it is running. If this happens, malfunctions can occur the next time the device is switched on.
- If you suspect a malfunction or faulty measured results, please inform the person responsible for the device immediately. If necessary, this person will contact the manufacturer or distributor to solve the problem.
- *Any serious incident that has occurred in relation to the device shall be reported to the manufacturer, the distributor and if applicable the competent authority of the European Member State in which the user and/or the patient is established.*
- *If there is any accidental skin contact with potentially infectious substances, e.g. human samples, disinfect the affected parts of the body with a suitable disinfectant solution.*
- *We expressly point out that if the equipment and consumables are not used correctly the protective measures provided for the equipment may be ineffective.*

4.3 Installing the device

Before start-up, check the supplied device and accessories for completeness referring to the list under 1.2.2. If anything is missing, please contact your supplier immediately.

Furthermore, please check all parts for intactness. Proper and safe operation is only guaranteed when using original parts and accessories. NEVER use damaged parts or parts from other manufacturers!

Place the device on a horizontal, smooth and dry surface *and within the permitted operating temperature range (see chapter 7 Technical data)*! Please choose a location where the device is protected from direct sunlight and extreme variations in temperature since these can impair measuring results.

Requirements on the set-up location

- no direct influence of humidity
- no direct sunlight
- no strong electromagnetic fields or ionizing radiation
- no rapid changes in temperature because of windows, doors, air-conditioning etc.
- level, water-proof support
- clearance over the entire footprint

Connecting the device to the power supply (see fig. 4.1 Power):

Please make sure that the voltage noted on the power adapter is the same as the voltage of your power grid.

The device is connected to the power supply via the included power supply adapter. Connect the power supply cable to the power supply adapter. Plug one end of the power supply adapter into the power connector at the right side of the device and the other end into the socket. **Note:** *The SUPER GL speedy is immediately in operation.*

Connecting EDP (see fig. 4.1 EDP):

Plug the EDP cable into the EDP-port at the right side of the device and connect the other end with your EDP. Please pay attention to the notes in the *EDP manual (available on request)* and also to the notes from your EDP partner.

The following image shows the interfaces on the right-hand panel of the casing of the SUPER GL speedy:

Fig. 4.1 Connections SUPER GL speedy (*EDP, Printer, Power*)

4.3.1 Installing accessory printer TPU-80

Overview of connections and operating elements on TPU-80

Fig. 4.2 TPU-80 printer

Operating elements (view from top):

- Side LED for offline indication (only shining through the housing, LED inside)
- LED
 - On: Printer is online
 - Off: Printer is offline (side LED on)
 - Flashes once briefly at intervals of 1 sec: out of paper (printer automatically offline, side LED on)
 - Flashes once for 1 sec at intervals of 1 sec: printer is offline, but data in buffer (side LED on)

 Feed → Paper feed

 On-/Offline → Start-/Stops printing

Fig. 4.3
Back side view TPU-80
(Interface; Mode; Power On/Off)

Connections (view from back side):

- Interface: For connection of printer cable
- Mode: Switch for setting different communications settings
- Power: For power supply connection
- On/Off: For switching on/off.

Connecting TPU-80 to SUPER GL speedy

If SUPER GL speedy is used together with the TPU-80 accessory printer, the printer is supplied with power via the SUPER GL speedy's integrated power connector, which is located on the back side.

The jack plug of the printer cable is inserted into the printer interface on the right-hand panel of the device's casing (see figure 4.1) and connected to the corresponding interface on the back panel of the printer (see figure 4.3).

When TPU-80 is turned on using the on/off switch for the first time, integrated short instructions will be additionally printed out when mode switch is in position 0.

In order to set the correct communication settings for SUPER GL speedy, please turn the mode switch to position 1. Selection becomes active after a restart of the printer.

Loading paper on TPU-80

To load the paper, pull the black lever on the cover forward and open the flap. In case of a new paper roll, remove the first layer, otherwise adhesive residues may contaminate / damage the printhead. Insert the paper so that it is guided from the bottom to the top at the front. Close the flap again and switch the printer to standby via the online key. The paper must always be torn off straight. Make sure that the tear-off edge is straight.

4.4 Initial operation of the device

When the device was installed as described above, *the following steps must be carried out for the device to start working:*

- Install the biosensor (chapter 6.3.2)
- Install the supply and waste containers (chapter 6.3.3)

Once this work has been completed, the installation of the device is finished. It runs through the necessary warm-up time and reaches the stand-by operating state.

To avoid loss of data, the device must not be switched off unless in “Stand by” mode or if required by a corresponding error message.

The switch on the right connection side is NOT an on/off switch!

The SUPER GL speedy is operated with the help of the Touch Screen only:

Fig. 4.4 Touch Screen

When operating the touch screen, please note:

- press only lightly
- do not use pointed or sharp objects
- do not use any solvents for cleaning, excepting the solution for cleaning and disinfection (chapter 6.3.1)

Browse menus and set numerical values by touching the respective areas. Please note, that dark shaped areas describe the status "OFF" or "INACTIVE", but light shaped areas describe the status "ON" or "ACTIVE".

While operating SUPER GL speedy you have to note, that there are some menu points / functions which you need for day-to-day operation and there are also some you only need in service cases.

The functions of day-to-day operation you could use without any special identification and they are changeable. The service function only should be used by special trained personnel because of it is protected by a keyword.

If there are not authorised contacts in protected areas the manufacturer assumes no liability for wrong measurement results!

The following page shows the menu structure of SUPER GL speedy. If you've got more questions, please ask the manufacturer or your distribution partner.


Fig. 4.5 Main program branches

4.5 Preparing the measuring process

4.5.1 Basics

The SUPER GL speedy uses prefilled reagents. For each analysis a prefilled reaction cup is needed.

For measuring with SUPER GL speedy you also have to use the *corresponding* biosensor, calibration *solution* and control material.

The cups you have to put into the *labelled* positions at the rotor and after that you can start to measure.

The following abbreviations have been used for marking the correct positions:

KA calibration solution

KO control solution

STAT quick sample

4.5.2 Sample preparation

Please observe the instructions on the package insert of the reagent cups concerning sample preparation!

The following notes complement the above notes and they are only valid if capillary blood is used as sample material:

When drawing a capillary blood sample, do not compress the tissue. Compressing the tissue leads to a dilution of the blood sample with intracellular fluid and can thus lead to faulty results. For taking a capillary blood sample, use suitable lancets and, if necessary, circulation-enhancing measures (such as lotions and massaging the spot) to yield a sufficient sample amount. *After sample taking, the capillary and the reaction cup may be contaminated with potentially infectious substances of human origin!*

On the following page, taking capillary blood using open-end capillaries is described and shown. Proceed likewise with an end-to-end capillary capillary (*this capillary has not to be broken*).



	<p>Taking capillary blood from the earlobe or the finger pad and filling the capillary over both markings</p>
	<p>Make sure, it is properly filled (sufficient amount of blood, no air bubbles, no drops of blood at the end of the capillary etc.)</p>
	<p>Carefully wipe off the outer surface of the capillary</p>
	<p>Break the capillary at the predetermined breaking point (predetermined breaking point is located in the middle between two markings)</p>
	<p>Insert the completely filled capillary into the pre-filled sample cup</p>
	<p>Shake sample cup until the blood has completely left the capillary</p>

Fig. 4.6 Sample preparation using open-end capillary

4.6 Measuring

4.6.1 Rotor operation

The SUPER GL speedy allows for measurements in two modes of operation.

- Determination of quick and control samples
- Determination of a sample series

To start a sample series proceed as follows:

- Place sample cups on the sample *rotor*. You do not have to start in position 1 and you may leave positions empty since the device possesses automatic sample cup recognition.
- In Stand by mode press START on the TOUCH SCREEN.

To measure a *quick* sample (STAT):

- Press STAT on the TOUCH SCREEN in Stand by mode or while measuring a sample series.

To start a control sample series proceed as follows:

- Only the occupied control positions are measured, no samples.
- Press the following areas on the TOUCH SCREEN in Stand by mode:
MENU → CONTROLS → MEASUREMENT

Prior to each of these measuring regimes, a calibration is carried out unless a different calibration option was chosen. Other calibration options are enabled and explained by service personnel upon customer request.

4.6.2 Chain operation

The determination of quick and control samples is carried out as described in 4.6.1.

To start a sample series proceed as follows (see p. 5.3.2.1):

- Switch on chain operation (Methods / Adjust / Count).
- Place sample cups on the sample *rotor*.

Chain operation is characterised by stopping the series at the first unoccupied position. You then have three ways to continue measuring samples:

- „Sample“ button
Measuring of the series is continued at the position flashing in the “Sample” area.
- „New rotor“ button
Measuring of the series is continued at position 1 on the *rotor*. The counting of the sequential numbers (numbers within the series) is continued with the next higher sequential number.
- „New series“ button
Measuring of the series is continued at position 1 on the *rotor*. Counting of the sequential numbers is restarted with number 1.

To measure a quick sample, press STAT on the TOUCH SCREEN in Stand by mode or while measuring a sample series.

To start a control sample series, i.e. only the occupied control positions are measured - no samples - successively press the following areas on the TOUCH SCREEN in Stand by mode:

MENUE → CONTROLS → *MEASUREMENT*

Chain operation is only recommended with an auto-calibration type. For more information see chapter 4.6.3.

4.6.3 Calibration

To guarantee correct measuring results the SUPER GL speedy needs a valid calibration. So the device calibrates for each *mode* of measurement automatically.

To customer request our service will activate also other calibration modes. In this case following areas are available by pressing menu item "Calibration":

Auto calibration (activated upon customer request):

Activates auto calibration mode. If the button is highlighted after touching, auto calibration is active.

Calibration takes place cyclically in the selected intervals.

In Stand by the period of time until the next calibration is displayed. In addition the "Cal off" button is displayed. Touching this button deactivates auto calibration.

Calibration stable (activated upon customer request):

This works just like auto calibration except that the next cyclic calibration is executed only when measuring is requested.

Start calibration

Executes a calibration.

If Auto calibration is activated, in the menu "Functions" - "Adjust" a new *button* "auto calibration" will be shown.

If the device is switched on and autocalibration mode is activated the device automatically starts the calibration after the warming time and remains in auto calibration mode.

Note: After switching on, the device will calibrate every time before the sample measurements for the first two hours.

4.6.4 Controls

The SUPER GL speedy has three control positions to meet the requirements of quality control.

To guarantee an effective quality control all three positions can be programmed.

Following input options are *currently available for the individual positions*:

- name of the control
- lower warning limit for glucose
- upper warning limit for glucose
- lower warning limit for lactate
- upper warning limit for lactate

By choosing the menu item "*Memory of controls*" the measured controls will be shown.

By pressing the button "Co1/Co2/Co3" the values for control 1, 2 and 3 *can be viewed or printed out selectively one after the other*. It will be printed the values from the dark shaped are up to the actual control value every time.

4.6.5 Method

The setting of the method depends on the used biosensor. The following settings can currently be made:

- Glucose
- Lactate
- Glucose and Lactate

4.7 Switching off the device

You may switch off the device, if there are no functions working. NEVER switch off the device during calibration or washing, otherwise the device may have malfunctions.

If the device is switched off for *more than 5 days or for transport see chapter 6.3.4 Decommissioning of the device / Disposal and 6.3.5 Prepare equipment for transport!*

Please store also the consumables (*especially sensors and controls*) as written at the package.

5 Operation - Part 2

5.1 Introduction

This part of the instruction manual describes special functions and settings relevant to the user. Furthermore, it gives additional information about quality control and some problems that can be solved by the user.

5.2 Menu functions

As described in chapter 4, there are two types of device functions: functions which are needed for day-to-day operation and functions which are only should be used by trained stuff.

For working with following functions you need besides special knowledge also exact knowledge about the menu structure of the SUPER GL speedy.

The overview of the menu structure you find in fig. 4.5.

5.3 Programming

5.3.1 Basics

The programming of controls, calibration *modes* and printer settings are described in chapter 4.

Following functions influence measuring results and their displaying and should only be activated by trained stuff (maybe by calling the service).

5.3.2 Method

5.3.2.1 Adjust

The following settings can be made here:

Cycle

- *defined the* measuring cycle
- Here you can select in which position of the series control measurements or recalibration take place.
- A cycle consists of a maximum of 8 segments.
- The „Back“ button deletes the last segment. The sample type is added by touching the button.
- In chain operation mode the cycle is executed until for the first time a sample is selected. Afterwards only samples are processed.

Sample repeat

Sample repetition - This function interacts with the “sample warning limits”. To select an item, touch the button. Highlighted buttons are active.

Count

- counting *samples*
- Here you can choose between a sequential number and a sample number of the day.
- Chain operation is intended for users needing devices with sample chains (see also 4.6.2).
- To switch chain operation on/off, select the menu item Chain operation. Two buttons (On, Off) are displayed. Touch the appropriate button. The active button is highlighted.

Reset parameters

This function resets the programmed parameters and the device returns into the delivery status. *All customer specific settings will be erased.*

In detail it happens as follows:

- sample number get number "1"
- time gets EU Standard
- year gets 4 digits
- language changed to German
- measuring unit gets mg/dl
- sample repetition gets status "OFF"
- sample number of the day gets status "OFF"
- calibration modus gets status "inactive" (device calibrates before each series and before each STAT)
- time of calibration stability gets 20 min
- time of auto calibration gets 30 min
- upper and lower warning limits of controls will be set to range limits
- sample warning limits will be set to range limits
- correction factor for methods gets standard values
- chain operation gets status "OFF"
- cycle gets standard value
- *EDP interface is set to standard*
- *barcode settings are set to standard*
- measuring value memory and control value memory were erased

Reset No. Of Samples

This function sets the actual sample number (number of the sequence or daily number) to number 1.

5.3.2.2 Select

This menu item *depends on the used sensor*. Here you can select the parameter which you want to measure, e.g. glucose or lactate as single parameter or both simultaneous from one sample.

5.3.2.3 Print

Prints out the parameter list.

5.3.2.4 Program

Unit

For the measuring units you can choose between mmol/l and mg/dl.

Warning limit

Enter a lower and an upper warning limit. For samples found outside these limits, an automatic repetition of the measurement can be programmed.

Value correction

Set measured value corrections

This function offers an *analyte separated* correction function with the equation $y = mx + n$ to adapt this device to a giant laboratory equipment, if it's necessary.

5.3.2.5 Plasmareference glucose

Using this function you are able to manipulate glucose values.

Using „Plasmaref. Glucose off“ the glucose values, which came from *whole* blood samples will be calculated as *whole* blood sample values.

Using „Plasmaref Glucose on“ the instrument will determinate the hematocrit values if the sample is a whole blood sample. Both values will be used to calculate the blood value into a plasma value. *If this function is activated, „Glucose-Hemolysate-Samples Plasmaref.“ appears in the printout* when the instrument starts a series or when it measures *quick* samples or controls.

5.3.3 Functions

5.3.3.1 Adjust

Adjust Time

Setting of the time in form HH:MM

Adjust Date

A changing of the date is not possible. *If necessary please contact our service hotline.*

Language

Currently are following dialog languages possible: German, English, Czech and Russian. More languages are possible after request. In this menu you also can choose the printing form of date and time.

EDP

Enables configuration of the EDP interface. This should only be done with *the knowledge of the EDP manual.*

5.3.3.2 Actions

Washing

Rinsing of the system

By the help of a pump solution will be pumped through the system from the suction hose in a defined time. This can be system solution or another solution, e.g. for disinfection. In this case the sensor will be destroyed.

Emptying

By the help of the pump the system will be emptied. In this case the suction hose should be outside the bottle with system solution.

Changing syringe

Used for maintenance (chapter 6.2)

Moves the plunger in a position for changing

Show datas of sensor

Shows the current data programmed in the sensor. In the display are shown following data:

```

sensor ok
samples:                xxxx
Remaining days          xx
method (example)       3 / 5 / 14 / 3 / 12
    
```

The digits mean in the sequence:

```

method number, optional reinforcing glucose, maximum reinforcing glucose, optional
reinforcing lactate, maximum reinforcing lactate
expirat. Time (Expiry)      month/year
lifetime                    months
activating                  date
    
```

Show supplies

In the display are shown following data:

```

supply ok
cups                      XXXX
buff.                    XXXX (buffer means system solution, specified in ml)
samples                  XXXXXX
Maintenance              DD.MM.YY
Sampl. to maint.        XXXXXXX (samples until next maintenance)
    
```

Adjust bubble detector

In the display 4 different values are shown.

By touching the key „Repeat“ an automatically incident for new adjusting the sensitivity starts.

After this incident the display shows *again 4 values* and an additional button “Save”. If there is a big difference between the data of the both upper lines, we recommend repeating this incident. With the button “Save” the actual values will be saved to adjust the bubble detector in the best way.

5.3.3.3 Read KEY CARD

Enables you to read the inserted Key Card. Supplies of system solution and sample cups are electronically stored on the Key Card. The Key Card thus guarantees the operation of the entire system (device and *consumables*) to insure the analysis quality.

A Key Card is included with every package of *System Solution Super GL and Glucocapil reaction cups*. By reading the Key Card, the supply of *system solution (buffer) or reaction cups (cups)* within the device is increased. When the supply is exhausted, the device reports "Please read KEY CARD". To read the Key Card, insert it into the *slot of the reader on the front side of the device*.

5.3.3.4 Service

For trained service personnel only. Description is in the service manual.

5.3.3.5 Barcode on/off

This function determines – separately for control, quick and normal samples *and calibration* – whether the bar codes on the cups are read or not. To switch, touch the button.

- Button dark = bar code is not read
- Button light = bar code is read

Note: The bar code reader is always active since it also serves to recognize occupied positions on the sample rotor.

If you want to use the bar code scanner to identify samples, please note the following:

- *The default* length of a barcode is 12 characters.
- Attach the barcode labels on the cups *at the correct height*. *When introducing a barcode label, check for any necessary adjustments!*
- The bar codes must be facing to the outside of the sample rotor.
- The device's operation is limited productivity-wise if sample cups without bar code labels are used while the bar code scanner is active.
- *For more information on the barcode, ask the service hotline for the EDP manual.*

Note: It is possible to improve the readability of particular bar codes by special configuration. In this case please contact the manufacturer or service.

5.3.3.6 User ID

This is used to enter a number for the operator using the device.

5.3.4 Measuring value memory

By choosing this menu item the measured rotor series will be shown.

The measuring value memory is a ring buffer. The most recent measuring values (approx. 500 values for one method; 250 values each by using two methods simultaneous) are stored. As soon as the buffer is full, the oldest values are deleted.

The measuring value memory *can only be completely erased* manually.

The stored values can be

- displayed
- sent repeatedly to the EDP
- printed repeatedly.

The Button "*calculate CV*" shows the mean value and the coefficient of variation (CV) of the last measured series.

To see the measuring value memory there are two display levels. In the level displayed first, all rows start with a „T“. A „T“ stands for a rotor, i.e. every row can contain a maximum of 60 samples. Use the arrow to select a rotor. If you press the „PRINT“ button, the results of the rotor marked in black are printed. Sending data to the EDP works correspondingly. If you press the „RESULTS“ button, the *single* results of the rotor marked in black are displayed. Now the values marked in black can be printed or sent.

5.3.5 Available cards

Key Card

Key Cards are supplied with packages of *System Solution SUPER GL and Glucocapil reaction cups*. They are different from each other and are used to read the supplies. Key Cards depend on the distributor.

6 Maintenance and troubleshooting

6.1 Introduction

This chapter gives information on maintenance of the SUPER GL speedy and about problems that may occur and how you might be able to solve them yourself.

If you are unsure about certain aspects, DO NOT try anything you might think appropriate without qualified technical help. DO NOT open the device without an authorised service technician! Please contact our service hotline!

6.2 Maintenance

The SUPER GL speedy needs maintenance once a year by *manufacturer-authorized and trained personnel*. *Authorised personnel are people who have gained expert knowledge by completing training courses offered by the manufacturer or authorised companies.*

The touch screen will print out a message after expiration of the maintenance rate.

Please contact the manufacturer or distributor immediately to make an appointment for service.

6.3 Servicing

The following operations can and should be carried out by the operator.

These actions are part of diligent care and serve to enhance the device's life span. They are NOT maintenance or service work, these may only be carried out by authorised service personnel!

6.3.1 Cleaning and disinfection

Please adhere to the regulations valid in your laboratory with regard to cleaning and disinfecting the device. For disinfection wipe the entire accessible surface of the device with a cloth containing disinfectant. Use a disinfectant for surface disinfection! Also note the instructions of the manufacturer of the disinfectant.

6.3.2 Exchange of sensor

Removal of Sensor

Open the sensor holder by pulling the lever

Fig. 6.1 Sensor holder closed

Remove the sensor

Fig. 6.2 Sensor holder open

Installation of Sensor

1. Open the packaging and take out the sensor
2. Open the sensor holder by pulling the lever (fig. 6.1)
3. Insert sensor
4. Close the sensor holder

6.3.3 Exchange of supply and waste container

To avoid interference with the operation of the device, it is recommended to exchange the containers only in Stand by mode and when the washing container and the sample canula are installed. The changing of these bottles should be done in a short time, because of washing the system there will be a little bit solution left from time to time. Furthermore, the supply container that was last emptied should be used as a waste container because the repeated use of containers could affect the function of the level sensors.

The waste connector and the waste container may be contaminated with potentially infectious substances of human origin!

Fig. 6.3 View of supply and waste container

Removal of supply and waste Container

1. Open the device cover
2. Remove the *connectors* (supply / waste hanging out on the back)
3. *Close the waste container*
4. Remove the containers

Insertion of supply and waste Container

1. If necessary, open the device cover
2. Insert the containers (supply *right*/ waste *left*)
3. Open the containers and insert the *connectors*
Keep the cap for closing the waste container later!
4. Close the device cover

6.3.4 Decommissioning of the device / Disposal

To decommission the device for a longer period (more than five days) or for transport, proceed as follows:

1. Rinse the device
2. Empty the system by removing the hose from the supply bottle and *repeat the rinsing procedure*.
3. Switch off the device and disconnect.
4. Please store all consumables as written at the packages.

Disposing of the device:

For the disposal of devices please ask your distributor.

6.3.5 Prepare device for transport

To prepare the analyser for transport, proceed as follows:

- *rinse and empty system*
- *switch off the device*
- *remove System Solution and waste container from analyser*
- *remove all cups from the analyser*
- *remove plug from power supply unit and printer and/or EDP*
- *Disinfect outside of the device.*
- *Pack and transport the analyser in the original cardboard box with foam insert, as this is the only way to avoid damage during transport as far as possible.*
 - *Place the device in one half of the foam, making sure that the analyser is placed correctly in the mould.*
 - *Place the other half of the foam on top and close the carton.*
 - *Attention, please! If force is needed to close the transport crate, the device is not correctly seated in the mould.*

A non conform preparation for transport can result in expensive damages at the device.

Please contact your distributor or the manufacturer if you have any further questions.

6.4 Error messages / Troubleshooting

6.4.1 Warnings

Before printing the measured results, the device checks whether set warning limits have been exceeded. The sample warning limits (chapter 3.2.3) are relevant for the determination of samples, control limits (chapter 4.6.4) are relevant for the determination of control samples. The following warning are displayed and printed:

Warning	Explanation
++++	Values above measuring range limit
-----	Values below measuring range limit
!!	Below or above the sample warning or control limits
*!	Previous control measurement outside the control limits and below or above the sample warning

Fig. 6.4 Table of displayed and printed warnings

6.4.2 Errors while measuring

Measuring value defective	Measuring signal is loaded with an interference signal which is more than 70 % of the measuring signal. The interference signal is compensated, but measuring value is less precise than a non disturbed signal. Correction value for plasma measuring is not plausible -> repeat measuring
Zero line <i>instable</i>	<ul style="list-style-type: none"> - appears incidentally -> repeat incident - changing reagent -> choose "Washing" - system obstructed -> rinsing by hand - sensor defective -> change sensor - electronic error -> call service
Zero comp instable	Analogue to error „zero line instable“
Calibration error <i>lower min. value</i>	<ul style="list-style-type: none"> - choosing wrong cup - sensor defective -> change sensor - electronic error -> call service
Calibration error drift to high	<ul style="list-style-type: none"> - appears incidentally -> calibrate new - rapid change of temperature (<i>draught</i>) -> calibrate new - system obstructed -> rinsing by hand - sensor defective -> change sensor - electronic error -> call service
Error sample	<ul style="list-style-type: none"> - system obstructed -> rinsing by hand - connection hose sample canula – sensor heavy discoloured - sample canula mal adjusted (defective)
Sample cup empty	- remaining quantity in the sample cup is not enough for a measurement

Fig. 6.5 Table of errors during the measuring

6.4.3 *Device errors*

If there are the following errors displayed the operator can't do anything without the help of the service personal. The displaying is only for the service, to describe more exactly the errors.

Type of errors

Communication error

ERROR (name of the component) no Ackn (no quittance of command)

Handler doesn't send a quittance of command

ERROR (name of the component) Timeout

Error of device / light barrier

ERROR (name of the component) light barrier doesn't open

Handler doesn't reach end position

ERROR (name of the component) 3 tries

ERROR (name of the component) 3 tries reference run

Problem of communication between controlling computers

ERROR (name of the component) device not available

Wrong parameter in the memory

ERROR (name of the component) Parameter-Flash

Hardware and Software don't fit together

ERROR (name of the component) Hardware version

Plunger doesn't reach the dependent end position

ERROR plunger type of error in clear text

Error lifter

ERROR lifter can't prick

Intern memory errors

ERROR Fe-RAM 1 (name of error)

There can be more similar messages of errors which are printed out to the same schema.

6.4.4 Measuring Errors

Analytes are determined using the enzymatic-amperometric measuring principle. The measuring signal is produced as a current change at an electrode due to a chemical reaction with the immobilized enzyme.

In the SUPER GL speedy a special measuring process is used that operates with a minimum amount of sample material. The hose between canula and sensor is critical for the proper functioning of this process. Therefore, this hose must only be replaced by an original spare part.

You can visually check proper operation by watching the air bubbles between sample and washing.

Sensor version *high speed*:

The air bubble that is drawn before a sample has to pass the sensor completely before sample will be measured.

The air bubble that is drawn after measuring the sample has to stop and stay for a while in front of the sensor.

Other sensor versions:

The air bubble that is drawn before a sample has to stop in front of the sensor until the sample taking arm has returned to the upward position.

The air bubble that is drawn after a sample has to stop in front of the sensor for a while.

As with all flow systems, unobstructed flow through and imperviousness of the canal between the sample canula and the hose pump is very important for the proper functioning of the device.

Always keep in mind that the hose pump works as a suction pump and that it produces a low pressure in the system. Leaking liquid is always a sign of permeability, e.g. worn seals, loose hoses or an incorrectly inserted sample canula.

Scattering measuring values for samples:

Could also be a result of incorrect sample preparation. Check by measuring a calibration cup several times. A defect sensor could also be the cause.

7 Technical data

Performance Data	
Measuring Time Series Measurement Glucose <i>high speed</i> Series Measurement Glucose / Lactate	ca. 18 seconds - 200 samples/hour ca. 30 seconds - 120 samples/hour
Measuring Range Glucose Lactate	0.5 mmol/l to 50 mmol/l or 9.01 mg/dl to 900.8 mg/dl 0.5 mmol/l to 30 mmol/l or 4.5 mg/dl to 270 mg/dl
Sample Volume	10 µl / 20 µl sample diluted with 500 µl / 1000 µl hemolysate solution
Precision (30 samples) - Variation Coefficient Indication Glucose (12 mmol/l or 216 mg/dl) Lactate (10 mmol/l or 90 mg/dl)	< 1.5 % < 1.5 %
Technical Specifications	
Working Temperature	+ 15 °C bis + 30 °C
Supply Voltage Power Consumption	24 V DC max. 60 W
<i>Classification according to Regulation 2017/746 on IVD medical devices (IVDR)</i>	<i>In vitro diagnostic medical device Risk class A</i>
Interfaces Printer EDP Barcode Reader	serial, RS 232 serial, RS 232 integrated
Dimensions / Weight <i>Width</i> <i>Hight</i> <i>Depth</i> Weight	470 mm 280 mm 450 mm 11 kg

Fig. 7.1 List of technical data