Troponin I control kit

for quality control of the Troponin I test kit parameter on the smart laboratory photometer.

EUROLyser

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Indication
Troponin I control kit

Kit Size
1 x 1ml

Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Troponin I control kit is to be used for internal quality control procedures in examining the accuracy and precision of quantitative Troponin I tests. For in vitro diagnostic use only! The controls are manufactured in a human serum base. The kit was tested by approved methods and found negative for the HIV I/II Ab, HBsAg and HCV. Because no method can offer complete assurance as to the absence of infectious agents, however, this material and all patient samples should be handled as though capable of transmitting infectious disease and disposed of accordingly.

Storage and Stability

Unopened vials are stable at 2 - 8°C until the expiration date on the label.

Once opened vials are stable until the expiry date when stored tightly capped at 2 - 8°C.

Preparation

Controls are liquid and ready for use. Take care to avoid the formation of bubbles. Use in the same way as a sample. Mix carefully before each use by swirling gently.

Dropper bottles:

Do not open the dropper bottle cap. Use 8 - 13 drops into an empty uncoated Eppendorf- or similar cup and pipette out of this cup.

Dispose the cup after use!

Round bottles:

Use only disposable pipette tips to avoid contamination.

Target Value and Range

Lot specific target value and range are stated in the value table. Results should be within the defined range. Troponin I concentration is expressed in ng/ml and µg/l.

Waste Management

Please refer to your local legal requirements.

Warnings and Precautions

Y Do NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbin to form explosive compounds.

Y Take the necessary precautions for the use of laboratory reagents.

Y Do not perform patient samples if the controls are outside accepted limits. If the result is not within the acceptable range repeat the measurement. If the result is again not within the acceptable range please contact your distributor.

References

1. Antman EM, Tanasijevic MJ, Thompson B, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996;335:1342-1349 | 2. Apple FS, Voss E, Lund L, Preese L, Berger CR, early detection of acute myocardial infarction and monitoring Chim Acta 1995;237:59-66 | 3. Mair J, Wagner I, Jakob G, Lechleitner P, Dienstl F, Puschendorfer B, Michel G: Different time courses of cardiac contractile proteins after acute myocardial infarction. Clin Chim Acta 1994;231:47-60 | 4. Mair J, Morandell D, Genser N, Lechleitner P, Dienstl F, Puschendorfer B: Equivalent early sensitivities of myoglobin, creatine kinase MB mass, creatine kinase isoform rations, and cardiac troponins I and T for acute myocardial infarction. Clin Chem 1995;41:1266-1272 | 5. Katus HA, Remppis A, Neuman FJ, Scheffold T, Diederich KW, Vinar G, Noe A, Matern G, Kuebler W: Diagnostic efficiency of troponin T measurements in acute myocardial infarction. Circulation 1991;83:902-912 | 6. Wang K, Asinger RW, Marriott HJ. ST-segment elevation in conditions other than acute myocardial infarction. N Engl J Med 2003;349:2128-2135

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
Referenzen

1. Antman EM, Tanasijevic MJ, Thompson B, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996;335:1342-1349 | 2. Apple FS, Voss E, Lund L, Preese L, Berger CR, early detection of acute myocardial infarction and monitoring Chim Acta 1995;237:59-66 | 3. Mair J, Wagner I, Jakob G, Lechleitner P, Dienstl F, Puschendorf B, Michel G: Different time courses of cardiac contractile proteins after acute myocardial infarction. Clin Chim Acta 1994;231:47-60 | 4. Mair J, Morandell D, Genser N, Lechleitner P, Dienstl F, Puschendorfer B: Equivalent early sensitivities of myoglobin, creatine kinase MB mass, creatine kinase isoform rations, and cardiac troponins I and T for acute myocardial infarction. Clin Chem 1995;41:1266-1272 | 5. Katus HA, Remppis A, Neuman FJ, Scheffold T, Diederich KW, Vinar G, Noe A, Matern G, Kuebler W: Diagnostic efficiency of troponin T measurements in acute myocardial infarction. Circulation 1991;83:902-912 | 6. Wang K, Asinger RW, Marriott HJ. ST-segment elevation in conditions other than acute myocardial infarction. N Engl J Med 2003;349:2128-2135

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
Processing of Troponin I control measurement

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






English

ATTENTION!

DO NOT wipe the pipetting tip after aspirating the sample! Allow single test min. 10 minutes to warm up to room temperature before use!

1. Preparation of test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press measurement button, enter required information using the touch screen

2. Select control

- 2.1 Sex menu: Select control
- 2.3 Instrument ready for start of analysis

3. Processing of Troponin I control measurement with pipette

- 3.1 Dispense some drops control liquid into serum cup. Aspirate 50µl from serum cup
- 3.2 Dispense 50µl sample INTO LIQUID in ERS cuvette
- 3.3 Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer
- 3.5 Start automatic sample processing by pressing the start button

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Technical details subject to change without notice.

Doc-ID: 874742b9-231a-46ad-b6ed-7b6409c9ff36 ; latest revision: DP 1.0 ; date of print: 2013-05-03