## BTA-Test (Blasentumorantigen) Test Komplementfaktor H (Immunchromatographie)

LD\_BG006-DE

20 Tests

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LD\_BG006

20 Tests

### FOR HEALTHCARE PROFESSIONAL USE ONLY

**PRODUCT NAME** LEPZI BTA Test (Bladder Tumour Antigen)

**SPEC:** LD\_BG006\_IFU\_001/Rev1.1

**Date of Revision:** 2023/09/28

**LEPZI BTA TEST (BLADDER TUMOUR ANTIGEN)  
(COMPLEMENT FACTOR H)**

LEPZI Bladder Tumour Antigen, BTA (Complement Factor H) Test is a point of care test for professional use which is an in vitro diagnostic test. This test is for **PROFESSIONAL USE ONLY** and allows users to perform tests using small sample volumes and to visually read the results quickly.

#### INTENDED USE

Bladder Tumour Antigen, BTA (also known as Complement Factor H) is a single use colloidal gold immunoassay intended for qualitative detection of BTA in urine. The test is intended as an aid in the assessment and diagnosis of patients on suspicion of urinary system lesions or bladder cancer in conjunction with cystoscopy and standard hospital procedures. The Lepzi BTA Test is intended for use by individuals trained in point of care settings. For professional use only and in patients  $\geq 18$  years of age.

#### SUMMARY AND EXPLANATION OF THE TEST

Cancer of the bladder, also known as urological cancer or urinary bladder cancer, is the 10th most common cancer in the world, and its incidence is steadily rising worldwide most common in the elderly, with smoking being a major risk factor (1). 90% of bladder cancer cases arise from these urothelial cells, mostly in the bladder but on rare occasions in the urinary tract (1). Bladder Tumour Antigen (BTA) in urine is a sensitive marker and has been used as an adjunct to cystoscopy in the diagnosis and follow-up of patients with bladder cancer (2).

#### PRINCIPLE

The BTA Test is an immunochromatographic assay utilising monoclonal antibodies to specifically detect the presence of BTA in urine. After successful sample addition, the BTA in the urine reacts with the colloidal gold-labeled mouse anti-human BTA antibody to form an immune complex. The immune complex migrates along the nitrocellulose membrane where it binds to fixed mouse anti-BTA antibody to form a visible line in the Test zone. If there is no BTA in the urine, no visible line will form in the test area. The sample continues to migrate to the control area where a visible line will indicate that the test has run correctly (Control line). The method has been standardised against a commercially available BTA Test.

#### COMPONENTS

BTA Test Components:

- The BTA Test strip contains mouse anti-human BTA monoclonal antibody (Test area) and goat anti-mouse IgG polyclonal antibody (Control area).
- The Test kit contains:
  - 20 BTA Test Cassettes
  - 20 disposable pipettes
  - 1 Instruction For Use (IFU) pack insert

#### MATERIALS REQUIRED BUT NOT PROVIDED

- Vortex Mixer
- Timer
- Specimen collection containers/cups

#### STORAGE AND STABILITY

- Store Test kit in a dry and cool (2~30°C) area out of direct sunlight.
- Shelf Life: 36 months. Do not use past expiry date printed on the pouch.

#### SPECIMEN

Urine.

#### STANDARDISATION

The Lepzi BTA Test has been standardised using control preparations of BTA referenced against a laboratory test.

#### WARNINGS AND PRECAUTIONS

- For use by healthcare professionals.
- For In-Vitro Diagnostic Use only.
- Tests should be performed between 15-30°C.
- Do not use after the expiry date of the test printed on the pouch.
- Keep the Test Cassette in the sealed pouch until ready to use and allow to come to room temperature. Use the Test Cassette immediately after removing from the pouch.
- Pipettes are for single use only. Discard after single use.
- Read at 5 minutes but no later than 6 minutes. Test result is not valid if read after 6 minutes.
- Discard and do not use any damaged or dropped Test Cassette or Pipette.
- Treat urine samples, Test Cassettes and used devices as if they were infectious.

#### TEST PROCEDURE

##### Attention:


- Urine samples should be collected into a clean, plastic urine cup (not paper or foam) free from preservatives.
- Specimens should be at room temperature prior to use if they have been stored refrigerated. Frozen or refrigerated urine samples must be returned to room temperature and thoroughly mixed before testing.
- Fresh urine samples should be gently inverted several times before testing to mix well.
- Urine samples should be tested as soon as possible after collection and within 4h when stored at 18-26°C. Samples can be kept refrigerated up to 48h. Frozen samples can be used up to 26 weeks.

##### Test Steps:

- Remove Test Cassette from pouch and place onto a flat, horizontal, clean surface.
- Label the Test Cassette with Patient Identification Number.
- Before testing, mix the urine sample thoroughly by inverting several times or using Vortex mixer.
- Take the Pipette and squeeze completely the top bulb, then insert the tip into the urine sample.
- Release the bulb slowly to fill the Pipette with the urine sample.
- Holding the Pipette over the sample well of the Test Cassette, add 3 full drops of urine without air bubbles to the sample well.
- Leave the test to develop for 5 minutes. Then read the result within 1 minute. See Interpretation of Results below.
- READ THE RESULT AT 5 MINUTES FROM SAMPLE APPLICATION BUT NO LATER THAN 6 MINUTES. AFTER 6 MINUTES FROM SAMPLE APPLICATION THE RESULT IS NOT VALID.**
- Discard Pipette and Test Cassette in biohazard container.

#### INTERPRETATION OF RESULTS

- Check the Control Area (C) position. A red or reddish-brown line must appear for the result to be valid.
- Positive Result:** A result is **POSITIVE** when a red or reddish-brown line is visible in both the Test Area (P) and the Control Area (C).
- Negative Result:** A result is **NEGATIVE** when a red or reddish-brown line is visible in the Control Area (C), but not in the Test Area (P).
- Invalid Result:** A result is **INVALID** when no red or reddish-brown line is visible in the Control Area (C).

#### PERFORMANCE CHARACTERISTICS

##### Clinical Sensitivity and Specificity

Lepzi BTA Test sensitivity and specificity was determined using urine samples from 245 patients collected from a hospital clinic. Both healthy patients with no symptoms and patients presenting with hematuria, and/or with symptoms of difficulty urinating and / or pain. All patients were assessed by cystoscopy. High risk patients with positive cystoscopy were also confirmed with biopsy to have urothelial cell pathology. Testing was completed at the clinic.

| LEPZI BTA Test | Clinical Diagnosis Results |          | Total |
|----------------|----------------------------|----------|-------|
|                | Positive                   | Negative |       |
| Positive       | 86                         | 32       | 118   |
| Negative       | 9                          | 118      | 127   |
| Total          | 95                         | 150      | 245   |

Sensitivity =  $(86/95) * 100\% = 90,53\%$

Specificity =  $(118/150) * 100\% = 78,67\%$

**Threshold Value:** The threshold value for a positive result is 22.4U/mL.

**Repeatability:** Using 10 Test Cassettes with a positive control gave 100% positive agreement.

**Inter-batch Difference:** 10 Test Cassettes were taken from each of 3 different batches of Test Cassettes. In total 30 Test Cassettes were tested with the positive control sample and demonstrated 100% agreement.

**Interfering Substances:** Samples containing up to 15mmol/L uric acid, 6.8 µmol/L bilirubin, 1g/L haemoglobin, and 10g/L albumin in urine samples were evaluated for potential cross-reactivity and interference. These substances did not impact the Complement Factor H (BTA) Test results.

#### LIMITATIONS

- Results should only be interpreted alongside cystoscopy and other standard of care procedures. This test is an aid in assessment and diagnosis only.
- The test result can be affected if the test is not brought to room temperature, if insufficient sample is added or if the sample contains non-specific interfering substances.
- The test is not validated for use in patients with kidney disease.

#### REFERENCES

- Saginala, K. et al (2020) Med.Sci.8, 15; doi:10.3390/medsci8010015
- Cheng, Z-Z. et al (2005) Clinical Chem 51: 856-863

#### CUSTOMER SERVICES

For product inquiries please contact Lepzi Customer Services at [customerservices@lepzi.com](mailto:customerservices@lepzi.com) or telephone **+(359)-897222462**. Any adverse results experienced with the use of this product, and/or quality problems should also be reported to Lepzi Customer Services on or by email: [customerservices@lepzi.com](mailto:customerservices@lepzi.com). For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In Vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

#### CONTACTS

Registrant/ manufacturer: Lepzi Diagnostika Ltd., Bulgaria  
Address: 111 "Tsarigradsko shose" blvd, Sofia, Bulgaria, 1784.  
Tel: +359-897-289-089  
After-sales Service: Lepzi Diagnostika Ltd., Bulgaria.  
Tel: +359-897-289-089  
Email: [customerservice@lepzi.com](mailto:customerservice@lepzi.com)

#### INTERPRETATION OF SYMBOLS


Consult instructions for use

For in vitro diagnostic use only

Do not reuse

Store at 2~30°C

CE Mark

Manufacturer

Sufficient for (n) tests

Reference Number

Not for Self-testing