# PT (INR) control kit

Control Kit for quality control for the PT (INR) parameter on the smart analyser.

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### Kit Size

2 x 1 mL

**Test kit preparation: Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

For use as accuracy control of the PT (INR) assay. Precision control, i.e. the recording of the day to day variation, is an important procedure in coagulation testing; checking the accuracy of the test system is particularly important. PT (INR) control kit is prepared from selected Austrian citrated plasmas ("Austrian Source") of at least 100 donations. The Controls contain stabilizers but no bactericide additives.

## Material Required

(not supplied with the kit)

- Pipettes
- Distilled water

## Storage and Handling

The expiry date printed on the labels applies to storage of the unopened bottles at +2 .. 8°C.

Stability after reconstitution:

- Room temperature: 4 hours
- Stored at +2 .. 8°C: 8 hours
- Stored at -20°C: 1 month

The vials can be only frozen once.

Upon storage, caps should be screwed tightly.

## Waste Management

Please refer to your local legal requirements.

## Test Procedure and Preparation

Open the vial carefully and reconstitute the contents in 1 mL of distilled water by carefully rotating the vial (avoid frothing).

Allow the reconstituted Coagulation Control N to stand for 10 minutes at room temperature before use.

Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

## Warnings and Precautions

IVD for in vitro diagnostic use

All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

Each single donor plasma and each lot of Coagulation Control N are tested and found negative for HbSAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

- Y Set the the samplotype to "Ctrl. Plasma" and sex to "control" before running the QC sample.
- Y Compare the results with the values indicated in the lot specific value sheet.
- Y If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor.

Control Kit zur Qualitätskontrolle des smart PT (INR)  
Test Kit am smart Photometer.

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## Processing of PT (INR) control measurement

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







## English

- 1. Preparation of test system**
  - 1.1 Insert RFID card
  - 1.2 Place ERS cuvette in smart test kit rack
  - 1.3 Place ERS cap in smart test kit rack
  - 1.4 Press symbol on touch screen, required information is entered by use of the touch display of the analyser
- 2. Attention: Select "Control" and "Ctrl. Plasma" – Measurement on the smart analyser**
  - 2.1 Sex menu: Select Control
  - 2.2 Sampletype menu: Select Ctrl. Plasma
  - 2.3 Instrument ready for start of analysis
- 3. Processing PT (INR) control measurement with pipette**
  - 3.1 Dispense some drops control liquid into serum cup. Aspirate 20 µl from serum cup.
  - 3.2 Dispense 20 µl sample into ERS cuvette - do not dispense below the line!
  - 3.3 Apply ERS cap firmly onto ERS cuvette
  - 3.4 Place ERS cartridge into smart analyser.
  - 3.5 Start automatic sample processing by pressing the .

### ATTENTION!

Do not wipe the capillary or pipetting tip after aspirating the sample! Allow single test min. 10 minutes to warm to room temperature!

## Wertetabelle für PT (INR) control kit Value sheet for PT (INR) control kit

g

### Control N

Target	value (min)	value (max)

H

### Control AK

Target	value (min)	value (max)

M Eurolyser Diagnostica GmbH  
www.eurolyser.com