hsCRP test kit

for quantitative in vitro determination of c-reactive protein in fingerblood, venous blood, serum or plasma on the smart or CUBE laboratory photometer.

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Kit size

32 tests/pack
1 x 1 ml (decision level)

Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

C-reactive Protein (CRP) is the best known among the acute-phase proteins, a group of proteins whose concentration increases in blood as a response to inflammatory disorders. A growing number of studies have examined that a more sensitive test, the so called hsCRP (high sensitive CRP) can predict recurrent cardiovascular disease, stroke and death in different settings. High levels of hsCRP consistently predict recurrent coronary events in patients with unstable angina and acute myocardial infarction (heart attack). Higher hsCRP levels also are associated with lower survival rates in these patients.

Method Measurement Range

Immunturbidimetric Assay:

Serum/Plasma: 0,25 - 20 mg/l

Whole blood: 0,50 - 20 mg/l

Principle

Kinetic determination of the concentration of CRP by photometric measurement at 546nm or 700nm of antigen-antibody reaction between antibodies to human CRP bound to polystyrene particles and CRP present in the sample.

Test Kit

ERS Cuvette pre-filled with buffer reagent
ERS Cap pre-filled with latex reagent with rabbit anti human CRP antibodies

Storage Instructions and Stability of the Test Kit

The reagents are stable until the indicated expiry date, if stored at 2-8°C.

Warnings and Precautions

The reagents contain sodium azide (0,95 g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Specimen Sample Material

Freshly drawn serum is preferred to guarantee highest test performance.
Heparin plasma, EDTA plasma, EDTA whole blood and capillary blood from finger tip may also be used if limited test performance is accepted.

Reference Range

Adults < 5mg/l (<0,5 mg/dl) ; values > 3mg/l (> 0,3 mg/dl) can indicate a risk for coronary heart disease CHD.
It is recommended that each laboratory establish its own expected range.

Quality Control

For internal quality control the hsCRP control kit should be used.
Order number: ST1005

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English

ATTENTION!

Do NOT wipe the pipette tip after aspirating the sample! Allow single test min. 10 minutes to warm up to room temperature before use!

1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press measurement button, enter required information using the touch screen

2. Select whole blood or serum

Please consult the laboratory photometer user manual for more information.

- 2.1 Sampletype menu: Select the sample type
- 2.2 Hematocrit correction (if desired)
- 2.3 Instrument ready for start of analysis

3. Processing for whole blood, serum or plasma

- 3.1a Aspirate 20µl blood from fingertip
OR ...
- 3.1b Aspirate 20µl venous blood from EDTA capillary blood collection system
OR ...
Aspirate 20µl venous blood from EDTA primary tube
OR ...
Aspirate 20µl serum from primary tube
OR ...
Aspirate 20µl plasma from centrifuged capillary blood collection system
- 3.2 Dispense 20µl sample INTO LIQUID in ERS cuvette
- 3.3 Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer
- 3.5 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

Technical details subject to change without notice.