D-Dimer control kit

Control kit for quality control for the D-Dimer parameter on the smart analyser.

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Kit size

16 Tests/Package
2 x 1 ml / high and low

Test kit preparation: Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

D-Dimer control kit is intended for use in the clinical laboratory as a control for quantitative determination of D-Dimer in the normal range and high range. Controls are prepared from human blood and each single plasma used for this lots are tested HBsAg, HIV 1/2 Ab and HCV Ab negative. Although handle the samples according to general rules in a clinical laboratory environment. Control kit is for for in vitro diagnostic (IVD) use only.

Stability and Storage

Store at 2-8°C. DO NOT FREEZE. Kit can be used until expiry date printed on the labels.

After reconstitution:

Stored at	Expiration
Room temperature	8 hours
15°C	8 hours
2°C	8 hours
-20°C	1 month

Warning

This testkit is for in vitro diagnostic use only. Do NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbin to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to your local legal requirements. All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

Preparation

The D-Dimer Controls are lyophilized and have to be reconstituted with 1 mL of 0,9% NaCl solution 10 minutes prior to use.

Note: Ones reconstituted there is no further preparation step necessary. (such as spinning or mixing with Citrat) It can directly be used same way like the patient sample!

Control values and Calibrator Standardisation:

No international calibrator is available for the standardization of D-Dimer. Therefore the D-Dimer calibrator is calibrated by using a commercial ELISA kit and the control values are assigned against the master calibration ELISA-based.

Accuracy should be checked according to your local/national QC requirements with D-Dimer Control high and D-Dimer Control low. The control plasma is treated in the same way as patient samples and results should be checked against the acceptable confidence limits supplied on the kit insert sheet and the vial label.

Warnings

- Y Measure the Controls same way like a patient sample and compare the values indicated in the lot specific value sheet included in the control kit.
- Y If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor.

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

Processing of D-Dimer control measurement

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



English

1. Preparation of test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette in smart test kit rack
- 1.3 Place ERS cap in smart test kit rack
- 1.4 Press symbol on touch screen, required information is entered by use of the touch display of the analyser

2. Settings on smart Photometer

- 2.1 In the gender/sex menu select "Control"
- 2.2 In the sampletype menu select "Citr. Plasma"
- 2.3 Instrument ready for start of analysis

3. Processing D-Dimer control measurement

- 3.1 Dispense some drops control liquid into serum cup. Aspirate 20µl from serum cup and pipette sample into ERS cuvette INTO LIQUID.
- 3.2 Apply ERS Cap firmly on ERS cuvette.
- 3.3 Place ERS cartridge into smart analyser
- 3.4 Start automatic sample processing by pressing the  start button.

ATTENTION!

Allow single test min. 10 minutes to warm to room temperature!

Wertetabelle für D-Dimer control kit Value sheet for D-Dimer control kit

g

Low Control

Target	value (min)	value (max)

H

High Control

Target	value (min)	value (max)

M

Eurolyser Diagnostica GmbH
www.eurolyser.com