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**SUPER GL** compact

Elektrochemischer Analysator

**REF** 0540 000 000 00

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Gültig ab 1. August 2022

## Hersteller / Copyright:

Dr. Müller Gerätebau GmbH  
Burgker Str. 133  
01705 Freital  
Germany

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  
**HITADO**

HITADO GmbH  
Dreihausen 2  
59519 Möhnese  
Germany

[www.hitado.de](http://www.hitado.de)

## Service – Hotline:

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[www.hitado.de](http://www.hitado.de)

## Service hotline:

If you have any questions concerning our products  
please contact us:

+49 (0)2924 9705-88

Fax: +49 (0)2924 9705-31

Email: [support@hitado.de](mailto:support@hitado.de)

**A) Table of contents:**

	Page
A ..... <i>Table of contents</i> .....	2
B ..... <i>Table of figures</i> .....	4
C ..... <i>Table of symbols</i> .....	5
1 ..... Preface .....	7
1.1 ..... Introduction .....	7
1.2 ..... The SUPER GL compact .....	7
1.2.1 ..... Basics .....	7
1.2.2 ..... <i>Device and accessories</i> .....	8
1.2.3 ..... Intended use .....	9
1.3 ..... Indication / Contraindication .....	9
1.4 ..... Manufacturers liability .....	10
2 ..... Safety .....	11
2.1 ..... Introduction .....	11
2.2 ..... Responsibility / Training of <i>the operator</i> .....	11
2.3 ..... General safety instructions .....	11
2.4 ..... Product-specific safety instructions .....	12
3 ..... Description of analyser .....	13
3.1 ..... Introduction .....	13
3.2 ..... Measuring principle .....	13
3.2.1 ..... <i>Electrochemical</i> measuring principle .....	13
3.2.2 ..... <i>Photometric</i> measuring principle .....	14
3.3 ..... <i>Calibration mathematics</i> .....	15
3.3.1 ..... <i>Calibration for electrochemical determinations</i> .....	15
3.3.2 ..... <i>Calibration for photometric determinations</i> .....	15
3.4 ..... <i>Consumables</i> .....	15
3.4.1 ..... <i>Containerkit SUPER GL</i> .....	16
4 ..... Operation - Part 1 .....	18
4.1 ..... Introduction .....	18
4.2 ..... Safety instructions .....	18
4.3 ..... Installing <i>the device</i> .....	19
4.3.1 ..... <i>Installing accessory printer TPU-80</i> .....	20
4.4 ..... <i>Initial operation</i> .....	21
4.5 ..... Preparing the measuring process .....	24
4.5.1 ..... <i>Basics</i> .....	24
4.5.2 ..... Sample preparation .....	24
4.6 ..... Measuring .....	26
4.6.1 ..... Calibration .....	26
4.6.2 ..... <i>Sample measurement</i> .....	26
4.6.3 ..... <i>Memory of results</i> .....	27
4.6.4 ..... Controls .....	28
4.6.5 ..... Method .....	28
4.6.6 ..... Plasma-related glucose values .....	29

	Page	
4.6.7	Function of the Operation Card	29
4.6.8	Operation with barcode reader	29
4.7	Switching off the device	30
5	Operation - Part 2	31
5.1	Introduction	31
5.2	Menu functions	31
5.3	<i>Configuration</i>	31
5.3.1	<i>Basics</i>	31
5.3.2	Method	31
5.3.2.1	<i>Choose</i>	31
5.3.2.2	Adjust	32
5.3.2.3	Program	32
5.3.2.4	<i>Show parameters</i>	32
5.3.3	Functions	32
5.3.3.1	Adjust	32
5.3.3.2	Actions	34
5.3.3.3	<i>Service</i>	34
5.4	<i>Control management</i>	35
5.4.1	General	35
5.4.2	<i>Control management in accordance with RiliBÄK</i>	35
5.4.3	<i>Programming of the quality control functions</i>	35
5.4.4	<i>Interpretation of the control memory</i>	37
6	Maintenance and troubleshooting	38
6.1	Introduction	38
6.2	<i>Maintenance</i>	38
6.3	<i>Servicing</i>	38
6.3.1	Cleaning and disinfection	38
6.3.2	<i>Exchange of the sensor</i>	39
6.3.3	<i>Exchange of the Containerkit</i>	39
6.3.4	<i>Decommissioning of the device / Disposal</i>	40
6.3.5	Prepare device for transport	40
6.4	<i>Error messages / Troubleshooting</i>	41
6.4.1	Warnings	41
6.4.2	<i>Errors while measuring</i>	41
6.4.3	<i>Device errors</i>	41
7	Technical data	42

**B) Table of figures**

	Page
<i>Overall view SUPER GL compact</i> .....	Fig. 1.1 ..... 7
<i>Package content</i> .....	Fig. 1.2 ..... 8
<i>Optional accessories</i> .....	Fig. 1.3 ..... 8
<i>Flow schema</i> .....	Fig. 3.1 ..... 13
<i>Reactions in the sensor for glucose</i> .....	Fig. 3.2 ..... 14
<i>Reactions in the sensor for lactate</i> .....	Fig. 3.3 ..... 14
<i>Schematic design of the photometer unit</i> .....	Fig. 3.4 ..... 14
<i>Backside of the device SUPER GL compact</i> .....	Fig. 4.1 ..... 19
<i>Connections SUPER GL compact</i> .....	Fig. 4.2 ..... 20
<i>TPU-80 Printer</i> .....	Fig. 4.3 ..... 20
<i>View of the backside TPU-80</i> .....	Fig. 4.4 ..... 21
<i>Touch Screen</i> .....	Fig. 4.5 ..... 22
<i>Main program structure</i> .....	Fig. 4.6 ..... 23
<i>Sample preparation using open-end capillary</i> .....	Fig. 4.7 ..... 25
<i>Sensor holder closed</i> .....	Fig. 6.1 ..... 39
<i>Sensor holder open</i> .....	Fig. 6.2 ..... 39
<i>Table of displayed and printed warnings</i> .....	Fig. 6.3 ..... 41
<i>Table of errors during measuring</i> .....	Fig. 6.4 ..... 41
<i>List of technical data</i> .....	Fig. 7.1 ..... 42

C) **Table of symbols**

**General used symbols**

Symbol	Description
	<i>Biohazard</i>
	Follow instruction <i>manual</i>
	In-vitro-Diagnostica
	Manufacturer
	<i>Distributor</i>
	<i>Date of manufacture</i>
	<i>For near-patient testing</i>
	CE Compliance
	Serial number
	Item number
	Contents of package
	Batch <i>number</i>
	Use <i>before</i>
	Storage temperature

**Symbols in instruction manual**


Symbol	Description
	<i>Attention</i>
<b>Bold</b>	<i>Very important notes</i>
italics	<i>Modifications</i>

Table of symbols

**Symbols on the display**

Symbol	Description
	Waiting period
	Blue reaction cup <i>for calibration</i> <i>prefilled with calibration solution for system glucose/lactate</i> <i>(if hemoglobin deactivated)</i>
	White reaction cup <i>for the sample</i> <i>prefilled with hemolysate system solution for system glucose/lactate</i> <i>(if hemoglobin deactivated)</i>
	Green reaction cup <i>for calibration</i> <i>prefilled with calibration solution</i> <i>for system glucose/lactate/hemoglobin</i> <i>(only if hemoglobin is activated)</i>
	Yellow reaction cup <i>for the sample</i> <i>prefilled with hemolysate system solution</i> <i>for system glucose/lactate/hemoglobin</i> <i>(only if hemoglobin activated)</i>
	Arrows upwards – remove reaction cup
	Arrows downwards – insert reaction cup
	Sample <i>cassette</i>
	Arrows to the right – remove sample <i>cassette</i>
	Arrows to the left – insert sample <i>cassette</i>
	Analysis time (calibration and sample <i>measurement</i> )
	Send to printer
	Send to EDP
	<i>Device locked due to invalid control measurement</i>
	Barcode reader connected
	<i>Measure single</i> sample with barcode reader
	<i>Measure cassette</i> with barcode reader

## 1 Preface

### 1.1 Introduction

*We congratulate you on your purchase of the SUPER GL compact analyser and wish you much success and pleasure in working with our analyser.*

*In the following chapter "The SUPER GL compact" you will get a first overview of your analyser: which analytes can be measured, which further devices and accessories belong to your analyser and an overview of the functioning of the device.*

*You will also receive information about safety instructions, liability and warranty, indications and contraindications of your analyser.*

*For more detailed and additional information, please read the further chapters.*

### 1.2 The SUPER GL compact

#### 1.2.1 Basics

The SUPER GL compact was *designed using the latest technology along with decades of experience in the area of production of clinical-chemical analysers.*

Fig. 1.1 Overall view SUPER GL compact

*It fulfils all legal specifications with regard to design and production that are required of all devices used in clinical laboratories. The compliance with the valid norms and statutes is documented through the visibly attached CE-Label. The CE-Label signifies compliance with all pertaining laws and regulations and consequently safety and confidence. Please refer to the corresponding Declaration of Conformity.*

*By employing an altogether newly developed technology for the determination of glucose and lactate, it is possible to fulfil all requirements of quality assurance (e.g. RiliBäk (Guidelines of the Federal General Medical Council for Quality Assurance in Medical Laboratories)) in medical laboratories while maintaining easy handling and minimum operating effort. All users are thus able to achieve analysis results that meet the quality demands.*

## 1.2.2 Device and accessories

### Deliverables:

<i>Designation</i>	<i>Quantity</i>
SUPER GL compact	1
Sample cassette	1
Power connection cable	1
Power supply adapter for device and printer	1
Container adapter set	1
Instruction manual	1

Fig. 1.2 Package content

### Optional

Accessory printer TPU-80 including printer cable	1
EDP cable (SUB-D9)	1
EDP cable (USB)	1
Barcode reader	1

Fig. 1.3 Optional accessories

### **1.2.3 Intended use**

*The SUPER GL compact is a fully automated analyser for professional biochemical in vitro diagnostics in clinical laboratories and laboratory-like environments for near-patient testing. Electrochemical biosensors, that work according to the enzymatic-amperometric principle can be used. In addition, photometric absorbance measurements can be performed at a fixed wavelength. The device must only be used and operated by trained personnel. Operating the device for self-testing is expressly forbidden.*

## **1.3 Indication / Contraindication**

### **Indication:**

The SUPER GL compact is an automatic analyser for the *determination of different biochemical analytes* in prediluted samples, e.g. in hemolysed blood samples. *Currently, biosensors are available for the quantitative measurement of glucose and lactate, as well as photometric determinations of hemoglobin. Therefore, the following chapters focus on these three analytes.*

*For the selection of sample material, please refer to package inserts of the biosensors and reaction cups. Use only the specified sample materials, else faulty measured results may occur.*

*The sample material is taken out of closed reaction cups which are placed on a sample cassette or sample rotor. The containers for the system solution and the waste solution are located as Containerkit at the backside of the device.*

*The device has the following features:*

- *automatic quick sample measurement from the integrated sample rotor or sample series measurement up to 6 samples from the sample cassette*
- *automatic calibration modes*
- *serial RS 232 printer interface*
- *serial RS 232 EDP interface*
- *PS2-interface for the external barcode reader*

### **Contraindication:**

Using *unsuitable* sample material can *result in faulty measuring results*. If in doubt, call the manufacturer!

Operating the device for *self-testing* is expressly forbidden!

## 1.4 Manufacturers liability

Legal liability and warranty claims are expressly excluded in the following cases:

- gross negligence or wilful damage of *the device*, parts *thereof* or *consumables*
- unauthorised opening of the *device* by untrained *personnel* (without *proper* service training)
- force majeure (e.g. *stroke* of lightning, water damage, fire)
- non-observance of instruction *manual* and package inserts

## 2 Safety

### 2.1 Introduction

The following chapters concern the safety of *the person operating* the device.

Read these chapters carefully **PRIOR** to *starting up the device* because they contain general safety warnings, warnings concerning the personal safety of the person operating the device, and warnings for the protection of the device.

***Displaying the following safety warnings does not release the person operating the device from adhering to the safety measures of the facility.***

### 2.2 Responsibility / Training of *the operator*

- *The SUPER GL compact analyser must only be used and operated by trained personnel. Trained personnel are people with an education with medical or laboratory scope and who have gained expert knowledge by studying this instruction manual and hands on at the SUPER GL compact. A training course by an employee of the manufacturer or of an authorised distributor can support gaining knowledge.*
- *Every user is responsible for adhering to safety, health and legal regulations, and operating the device only according to its intended use.*
- *Interpreting the results and diagnosing on that basis must be left to a medical practitioner. Operating the device for self-testing is expressly forbidden.*
- *The operator must ensure the continuous validation of the correct function of the device through effective quality assurance measures. To this end, the local and national quality assurance requirements must be observed. At least after installation, after each new start-up and after troubleshooting, control measurements must be carried out with a suitable control material and checked for correctness.*

### 2.3 General safety instructions

- *Prior to operating the device, read the entire instruction manual especially the instructions for sample taking. If you have any questions, please contact the manufacturer or authorised distributor.*
- *Every person working with the device must be acquainted with the relevant safety rules prior to operating the device and these rules should be kept at hand all the time.*
- *Please pay attention to all general safety rules for laboratories such as wearing protective gloves, and the applicable disinfection and hygiene regulations.*
- *If there is any accidental skin contact with potentially infectious substances, e.g. human samples, disinfect the affected parts of the body with a suitable disinfectant solution.*
- *To avoid risk of electric shock, do not place the device or power supply in water or other liquids! If the cable or the power supply adapter is damaged in any way, you must not continue using the power supply. Never touch the plug of the power supply adapter with wet hands. The power supply adapter must only be used indoors and must be protected from humidity.*

## 2.4 Product-specific safety instructions

- *The device may only be used for the intended use with special attention to the defined usage restrictions and constraints that have to be strictly adhered to (consult the manufacturer if necessary).*
- *Operate the device only on smooth, horizontal surfaces. Avoid variations in temperature, drafts, direct sun light, and vibrations. These can result in faulty measuring values.*
- *In case of malfunctions, stop operating the device immediately! Prior to continuing to operate the device, read the notes concerning cleaning, error messages and troubleshooting. After consulting the manufacturer or authorised distributor you may ship the device for repairs to the manufacturer or authorised distributor.*
- *Use only original accessories and spare parts to avoid damage to device and people. Repairs must only be conducted by the manufacturer or by companies authorised by the manufacturer!*
- *The use of reagents and consumables that are not expressly recommended by the manufacturer can cause severe measuring errors and malfunctions and is therefore not permissible.*
- *If the user opens the device without authorisation, the user shall not be entitled to any rights concerning the liability for the device and damages caused thereby.*

### 3 Description of the analyser

#### 3.1 Introduction

This chapter describes the *analyser's* measuring principle, *the calibration mathematics and the consumables*. For detailed instructions and descriptions of the device, please refer to chapter *Operation*.

#### 3.2 Measuring principle

##### 3.2.1 Electrochemical measuring principle

*The determination of analytes* with the SUPER GL compact is based on an electrochemical measuring principle with a biosensor. *With the help of* an analyser pump system solution, calibration solution, control material or patient's material is conveyed through a sensor. The electrodes inside the sensor are separated from the liquid by semipermeable layers containing *e.g.* immobilized enzymes. The following *figure* shows the flow *schema*. *The electrochemical reaction takes place inside the biosensor*.


Fig. 3.1 Flow *schema*

*The next two figures* show these reactions as example for the measurement of glucose and lactate.


Fig. 3.2 Reactions in the sensor for glucose


Fig. 3.3 Reactions in the sensor for lactate

### 3.2.2 Photometric measuring principle

The photometer unit, consisting of LED, flow cuvette, interference filter, photo detector and electronic evaluation unit, is arranged in front of the biosensor and the pump. These are used to determine an extinction equivalent value for the sample solution. The dependence of the extinction on the concentration is described by the Lambert-Beer-Bouguer law. Under the constraints defined there, the size of the extinction is proportional to the material concentration.


Fig. 3.4 Schematic design of the photometer unit

### **3.3 Calibration mathematics**

#### **3.3.1 Calibration for electrochemical determinations**

The change of the electrical current on the biosensor will be used to calculate the concentration of the analyte based on a calibration curve. Because the correlation between concentration of analyte and measuring signal is linear a linear calibration mathematics is used. The concentration of the sample will be calculated by the following equation:

$$C_{Sample} = \frac{C_{Calibrator}(I_{Sample} - I_{SL})}{I_{Calibrator} - I_{SL}}$$

$C_{Sample}$  = Concentration of sample

$C_{Calibrator}$  = Concentration of calibrator

$I_{Sample}$  = Current of sample

$I_{Calibrator}$  = Current of calibrator

$I_{SL}$  = Current of System Solution

#### **3.3.2 Calibration for photometric determinations**

The photometric absorbance measurement is used to calculate the concentration of the analyte using an empirically determined extinction factor. Since the correlation between the concentration of the analyte and the measured absorbance is linear, linear calibration mathematics is used. The concentration of the sample is calculated using the following equation:

$$C_{Sample} = F_{EX}(A_{Sample} - A_{Blank})$$

$C_{Sample}$  = Concentration of sample

$F_{EX}$  = Factor of extinction

$A_{Sample}$  = Absorption of sample

$A_{Blank}$  = Absorption of system solution

### **3.4 Consumables**

To operate the analyser following consumables are needed:

- Glucocapil reaction cups prefilled with Hemolysate-System-Solution including Operation Card
- capillaries or pipette for dosing sample
- Glucocapil calibration solution
- Containerkit SUPER GL with 300ml System Solution SUPER GL
- biosensor for the analyte
- control material

For detailed instructions on the use of these consumables, please refer to the following chapters of this manual respectively the separate product information of these consumables.

### 3.4.1 Containerkit SUPER GL

A combination of system solution / waste container made specifically for the SUPER GL compact can be attached to the back of the device. The sizes are such that the waste container is full when the system solution is used up. The device monitors the system solution but not the waste container.

The System Solution SUPER GL is intended for use as washing solution for automatic filling and washing the liquid system of the analyser SUPER GL compact and to create a stable working environment for the biosensors.

The System Solution SUPER GL is not corrosive, poisonous or seriously harmful to skin. If there is any skin contact with the System Solution, rinsing with water is sufficient.

The System Solution SUPER GL must be stored protected from light and at temperatures of +2°C and above. It is ready to use and stable until the indicated expiry if protected from contamination, proper storage and handling are observed.

#### Installation Containerkit

1. Make sure that the cap of the filled container is white. The use of Containerkit with a different lid colour can lead to incorrect measurement results!

2. Hang the Containerkit on the hooks on the back of the device and unscrew the cap of the filled container. Keep the cap for closing the waste container later!


3. Screw the white connector onto the filled container and the red connector onto the empty waste container!

### Changing Containerkit

*The Containerkit may only be exchanged as a whole. Never fill up the supply container!*

1. Remove both connectors from Containerkit!

***The red connection and the filled waste container may be contaminated with potentially infectious substances of human origin!***

2. Close the filled waste container by using the white cap which was kept by installation.

3. Remove the Containerkit from the hooks on the back of the device!
4. Install a new Containerkit as described above!

***Do not use other than the prescribed bottles or fill up the system solution container since this can result in malfunctions.***

*For disposal, system solution and waste containers can be separated from each other by breaking the label. Please refer to national and local legal requirements.*

## **4 Operation - Part 1**

### **4.1 Introduction**

In this part of the instruction *manual*, all information is *included that is useful for the day-to-day operation of the device*.

In part 2, all additional information is *included that is important for understanding the functions, complementary functions and certain sources of possible problems*.

*The qualified personnel for the device must be familiar with both parts and must also have the medical knowledge to be able to interpret the acquired values correctly. Drawing conclusions for a diagnosis and a therapy is reserved exclusively for the medical specialist.*

### **4.2 Safety instructions**

As mentioned before, *certain safety warnings must be heeded when operating the device to guarantee correct and faultless operation*:

- *The device must only be used for the described intended use and indication and must only be used and operated by trained personnel.*
- *Every user is responsible for adhering to safety, health and legal regulations, and operating the device only according to its intended use.*
- *Interpreting the results and diagnosing on that basis must be left to a medical practitioner. Operating the device for self-testing is expressly forbidden.*
- *In daily operation, regular checks of the results should be made; if needed, an additional control measurement should be carried out.*
- *Do not switch off the device or disconnect it from the power supply while it is running. If this happens, malfunctions can occur the next time the device is switched on.*
- *If you suspect a malfunction or faulty measured results, please inform the person responsible for the device immediately. If necessary, this person will contact the manufacturer or distributor to solve the problem.*
- *Any serious incident that has occurred in relation to the device shall be reported to the manufacturer, the distributor and if applicable the competent authority of the European Member State in which the user and/or the patient is established.*
- *If there is any accidental skin contact with potentially infectious substances, e.g. human samples, disinfect the affected parts of the body with a suitable disinfectant solution.*
- *We expressly point out that if the equipment and consumables are not used correctly the protective measures provided for the equipment may be ineffective.*

### 4.3 *Installing the device*

*Before start-up, check the supplied device and accessories for completeness referring to the list under 1.2.2. If anything is missing, please contact your supplier immediately.*

Furthermore, please check all parts for intactness. Proper *and safe* operation is only *guaranteed when using* original parts and accessories. NEVER use damaged parts or parts from *other manufacturers!*


*Place the device on a horizontal, smooth and dry surface and within the permitted operating temperature range (see chapter 7 Technical data)! Please choose a location where the device is protected from direct sunlight and extreme variations in temperature since these can impair measuring results.*

#### Requirements on the set-up location

- no direct influence of *humidity*
- no direct sunlight
- no strong electromagnetic fields or ionising radiation
- no rapid changes in temperature because of windows, doors, air-conditioning etc.
- level, water-proof support
- *full ground clearance required over the entire installation area*

#### Installation container adapter set:

*Connect the hose of the red adapter to the hose connection of the hose marked with a red ring on the back of the unit! Connect the hose of the white adapter to the hose connection of the second hose on the back of the unit! Fix both hoses in the hose holder!*

*Fig. 4.1 Backside of the device SUPER GL compact*

#### Connecting the device to the power supply (s. Fig. 4.2):

Please make sure that the voltage noted on the power *adapter* is the same *as the voltage of your power grid.*

*The device is connected to the power supply via the included power supply adapter. Connect the power supply cable to the power supply adapter. Plug one end of the power supply adapter into the power connector at the right side of the device and the other end into the power socket.*

**Connecting EDP (s. Fig. 4.2):**

Plug the EDP cable into the EDP-port at the right side of the device and connect the other end with your EDP. Please pay attention to the notes in the EDP manual (available on request) and to the notes from your EDP partner.

The following image shows the interfaces on the right-hand panel of the casing of the SUPER GL compact:


Fig. 4.2 Connections SUPER GL compact

**4.3.1 Installing accessory printer TPU-80**

Overview of connections and operating elements on TPU-80

Operating elements (view from top):

- Side LED for offline indication (only shining through the housing, LED inside)
- LED
  - On: Printer is online
  - Off: Printer is offline (side LED on)
  - Flashes once briefly at intervals of 1 sec: out of paper (printer automatically offline, side LED on)
  - Flashes once for 1 sec at intervals of 1 sec: printer is offline, but data in buffer (side LED on)

➤➤ Feed → Paper feed

⏸ On-/Offline → Start-/Stops printing

Fig. 4.3 TPU-80 printer

**Fig. 4.4**  
*Back side view TPU-80  
(Interface; Mode; Power On/Off)*

*Connections (view from back side):*

- Interface: For connection of printer cable*
- Mode: Switch for setting different communications settings*
- Power: For power supply connection*
- On/Off: For switching on/off*

#### Connecting TPU-80 to SUPER GL compact

*If SUPER GL compact is used together with the TPU-80 accessory printer, the printer is supplied with power via the second connection of the power supply adapter. If necessary, remove the cap from the second plug of the power supply adapter.*

*The jack plug of the printer cable is inserted into the printer interface on the right-hand panel of the device casing (see figure 4.2) and connected to the corresponding interface on the back panel of the printer (see figure 4.4).*

*When TPU-80 is turned on using the on/off switch for the first time, integrated short instructions will be additionally printed out when mode switch is in position 0.*

*In order to set the correct communication settings for SUPER GL compact, please turn the mode switch to position 4. Selection becomes active after a restart of the printer.*

#### Loading paper on TPU-80

*To load the paper, pull the black lever on the cover forward and open the flap. In case of a new paper roll, remove the first layer, otherwise adhesive residues may contaminate / damage the printhead. Insert the paper so that it is guided from the bottom to the top at the front. Close the flap again and switch the printer to standby via the online key. The paper must always be torn off straight. Make sure that the tear-off edge is straight.*

## **4.4 Initial operation**

*When the device was installed as described above, proceed as follows:*

- Install the biosensor (chapter 6.3.2)*
- Install the Containerkit (chapter 3.4.1)*
- Install the Operation Card (chapter 4.6.7)*

Once these tasks have been completed the installation of the equipment is finished. The equipment is switched on by activating the switch.

*To avoid loss of data, the device must not be switched off unless in “Standby” mode or if required by a corresponding error message.*

After the equipment has been switched on, it will run in for the necessary time and will then be ready for work.

*The SUPER GL compact is operated exclusively by means of the touch screen:*

*Fig. 4.5 Touch Screen*

When *operating* the touch screen, please note:

- *press only lightly*
- *do not use pointed or sharp objects*
- *do not use any solvents for cleaning, excepting the solution for cleaning and disinfection (chapter 6.3.1)*

Leafing through the menu or adjusting the numbers is carried out by touching the relevant buttons. Note that dark underlined buttons describe the “off” or “inactive” state, while light underlined or blinking buttons show the “on” state.

When operating the SUPER GL compact note that there are menu points/functions which may be needed for daily work and those which are only used for servicing.

Functions which are needed on a daily basis can be called up without special passwords and can also be changed. However, the service functions are protected with a password known only to authorised staff who are the only ones allowed to use it.

***If there are not authorised contacts in protected areas the manufacturer assumes no liability for wrong measurement results!***

*The following page shows the menu structure of SUPER GL compact. If you’ve got more questions, please ask the manufacturer or your distribution partner.*


Fig. 4.6 Main program structure

## 4.5 Preparation of measuring process

### 4.5.1 Basics

The SUPER GL compact uses *prefilled* reagents. *For each analysis a prefilled reaction cup is needed. For measuring with SUPER GL compact you also have to use the corresponding biosensor, calibration solution and corresponding control material.*

The cups are placed in the *correspondingly marked* positions on the sample rotor or the sample cassette and after that you can start to measure.

*The following colour fields were used to mark the correct positions on the sample rotor:*

*Blue calibration solution*

*Red single sample*

*Coloured rings (orange and red) control position*

*(When using GL Control N and GL Control P, the colour of the cups corresponds to the colour rings!)*

Measuring starts automatically if a sample cup is placed in the sample position (red segment) on the sample rotor. Calibration is carried out if necessary, before measurement.

**Note:** **Cassette must not be positioned in the analyser.**

Series measurement starts automatically if the cassette is positioned in the analyser. Only the occupied positions will be processed.

A control measurement must be indicated on the touch screen. On instruction, the rotor is turned to the loading position for the control samples. After that, the positions must be loaded within 15 seconds and the measurement started by instruction *again* by using the touch screen. If this does not happen, the analyser reverts to the initial state.

### 4.5.2 Sample preparation

Please observe the instructions on the package insert of the reagent cups concerning sample preparation!

*The following notes complement the above notes and they are only valid if capillary blood is used as sample material:*

***When drawing a capillary blood sample, do not compress the tissue. Compressing the tissue leads to a dilution of the blood sample with intracellular fluid and can thus lead to faulty results. For taking a capillary blood sample, use suitable lancets and, if necessary, circulation-enhancing measures (such as lotions and massaging the spot) to yield a sufficient sample amount. After sample taking, the capillary and the reaction cup may be contaminated with potentially infectious substances of human origin!***

On this page, taking capillary blood using open-end capillaries is described and shown. Proceed *likewise* with an end-to-end capillary (*this capillary has not to be broken*).



	<p>Taking capillary blood from the earlobe or the finger pad and filling the capillary over both markings.</p>
	<p>Make sure, it is properly filled (sufficient <i>amount</i> of blood, no air bubbles, no drops of blood at the end of the capillary etc.)</p>
	<p>Carefully wipe off the <i>outer surface</i> of the capillary <i>Make sure that the capillary is still filled over both markings!</i></p>
	<p>Break the capillary at the <i>predetermined breaking point</i> (<i>predetermined breaking point is located in the middle between two markings</i>)</p>
	<p>Insert the completely filled capillary into the prefilled <i>reaction cup</i></p>
	<p>Shake the reaction cup strongly until the blood has completely dissolved from the capillary.</p>

Fig. 4.7 Sample preparation using open-end capillary

## 4.6 Measuring

### 4.6.1 Calibration

To ensure correct *measurement* results, the SUPER GL compact requires a valid calibration. *The SUPER GL compact performs the calibrations independently and monitors their validity. A calibration is valid for a maximum of 60 minutes. The validity period may vary due to usage and environmental conditions.*

*For the calibration, you can choose between two calibration variants. An explanation can be found in chapter 5.3.3.1 under "Autocal".*

*Before the first calibration and after an unsuccessful calibration, the device requires the insertion of a new calibration cup. The calibration cup must always be placed on the sample rotor in the blue segment.*

***The colour of the calibration cup to be used in the display must be observed.***

***Even if the calibration solution in a cup is not used up in one working day, a new calibration cup must be used every day.***

*The start of a calibration can also be initiated manually in the device menu under Functions/Calibration.*

### 4.6.2 Sample measurement

With the SUPER GL compact it is possible to carry out measurements in two types of operation.

- determination of single or control samples
- determination of a sample series

*The following applies to all operating modes:*

If a valid calibration is not available, a blue or green reaction cup is displayed on the touch screen with an arrow (empty position: downward arrow; loaded position: upward arrow). After insertion of a calibration cup the analyser carries out a calibration.

After calibration has been carried out or if a valid calibration is available, a reaction cup for a blood sample and a cassette symbol with appropriately positioned arrows and the last measured result are displayed on the touch screen.

***The colour of the reaction cups to be used in the display must be strictly observed!***

To determine a single sample, proceed as follows:

- *Insert* a sample into the red segment on the *sample rotor*. *The measurement starts automatically.*

To start a control sample series, proceed as follows:

- For a control sample series only the occupied control positions are measured and no *other* samples.
- A control *measurement* must be indicated on the touch screen. On instruction, the *rotor* is turned to the loading position for the control samples. After that, the positions must be loaded within 15 seconds. *The orange ring refers to control 1 set in the device and the red ring to control 2. (When using GL Control N and GL Control P, the colour of the cups corresponds to the colour rings!)*
- *The measurement must be started by instructing again* using the touch screen. If this does not happen, the analyser reverts to the initial state.

To start a sample series, proceed as follows:

- Load the *cassette* with reaction cups. It is not necessary to begin with position 1 and positions can be left empty, as the analyser has automatic reaction cup recognition.
- Place the *cassette* in the analyser. Measuring begins automatically. During series measuring it is not possible to scan *single* samples or control samples.

### **4.6.3 Memory of results**

The data memory is designed as a ring memory. It stores the results of the last *600 measurements (single, control, series measurement)*. If it is full, the oldest value is deleted.

The stored values can be

- *viewed*
- repeatedly sent to the EDP
- repeatedly printed
- deleted.

The data memory is permanently available even though no *measurement* is being carried out. The last result measured is displayed. There is an arrow on the right which indicates upwards. By clicking on this arrow the next oldest result is displayed. Two arrows are now displayed, up and down.

*By touching the printer symbol at the bottom right, either the measured value memory from the displayed sample to the last measured value or the displayed individual value is printed out again, depending on the selection. The same applies to the EDP symbol, whereby only the currently displayed data record can be sent to the EDP.*

*The data memory can also be selected manually in the device "Menu" under "Memory of results". Here you have the option of completely deleting the measured value memory manually.*

#### 4.6.4 Controls

The SUPER GL compact has two control positions to meet quality assurance requirements.

To ensure effective quality assurance, the two control positions can be individually programmed. Control 1 is the position with the orange-coloured ring and control 2 is the position with the red ring. The following input options are *currently available* for the individual positions:

- name of control
- lower warning limit for glucose
- upper warning limit for glucose
- Lower warning limit for lactate
- Upper warning limit for lactate
- Lower warning limit for haemoglobin
- Upper warning limit for haemoglobin

For the control positions, the stored values can always be displayed *in the memory of results*. *The separate control memory is only active when using the quality control function.*

You can find further information on control *management* in chapter 5.4.

#### 4.6.5 Method

*The choice of methods depends on the biosensor used. The following are currently available:*

- Glucose / Glucose plasma-related
- Lactate
- Haemoglobin

*A combination of the methods is possible.*

**For the measurement of hemoglobin special calibration cups (green) and reaction cups (yellow) must be used.**

#### **4.6.6 Plasma-related glucose values**

*Leading professional organisations for clinical chemistry and diabetes recommend that results of glucose measurements be reported exclusively as glucose in plasma, regardless of which sample material was used for glucose measurement.*

*By activating the function "GL-Plasma-rel on" it is possible with the SUPER GL compact to output glucose measurements performed with whole blood samples as glucose in plasma.*

*If this function is activated, the haematocrit value is also determined internally in addition to the glucose. With the help of the determined haematocrit value, an individual conversion factor is calculated for each sample and the measurement result is output plasma-related as glucose in plasma. The result for plasma-related glucose measurement values from whole blood samples can be recognised by the measurement value identification GLU-P.*

*If sample materials other than blood are used when the function "GL-Plasma-rel on" is activated, the measurement result is not converted because the device recognises this by the internal determination of the haematocrit value. Results that are not converted to glucose in plasma can be recognised by the measurement value label GLU.*

*This ensures that when the "GL-Plasma-rel on" function is activated, all results are output as glucose in plasma, regardless of whether whole blood or plasma is used as the sample material.*

#### **4.6.7 Function of operation card**

Stocks of reaction cups are electronically stored on the operation card. The operation card therefore guarantees the use of the full system (device and consumable material) to ensure analysis quality. *Operation Cards depend on the distributor.*

The relevant operation card is enclosed with each package of *Glucocapil prefilled* reaction cups. During operation, the Operation Card must be inserted *into the card slot on the top of the SUPER GL compact with the imprint facing forwards.*

#### **4.6.8 Operation with barcode reader**

An external barcode reader can be connected to the SUPER GL compact. The machine recognises when a barcode reader is connected and waits for a read process for each sample before measurement. The barcodes are assigned to the relevant sample positions during measurement.

If you want to use the barcode reader for sample identification, note the following instructions:

- The maximum length of a barcode is *13 characters*.
- *Under Menue/Length Barcode, the length of the barcode used for data transmission and display can be defined. This allows you to increase the security for your unique sample identification.*

When operating the *device* the symbols (reaction cup and sample *cassette*) are appropriately labelled on the screen. Note that a red margin around the appropriate symbol indicates that a barcode has not yet been read and the *device* can therefore not measure the sample. Conversely, a green margin around the relevant symbol means that at least one barcode has been measured and the machine is ready to measure the relevant sample or sample *cassette*.

#### **4.7 Switching off *the device***

*You may switch off the device, if there are no functions working. NEVER switch off the device during measurement, calibration or washing, otherwise the device may have malfunctions.*

*If the device is to remain switched off for more than 5 days or is to be prepared for transport, follow chapters 6.3.4 Decommissioning of the device / Disposal and 6.3.5 Preparing the device for transport.*

*Please store also the consumables (especially sensors and controls) as written at the package.*

## **5 Operation - Part 2**

### **5.1 Introduction**

This part of the instruction manual describes special functions and settings relevant to the user. *Furthermore, it gives additional information about quality control and some problems that can be solved by the user.*

### **5.2 Menue functions**

As described in chapter 4, there are two types of device functions: functions which are needed for *day-to-day* operation and functions which are only should be used by *trained* stuff.

For *working with* following functions you need *besides special* knowledge also exact knowledge about the menu *structure* of the SUPER GL compact.

The overview of the menu structure you find in fig. 4.6.

### **5.3 Configuration**

#### **5.3.1 Basics**

*On SUPER GL compact you have the option of making customer-specific settings with regard to the measuring operation already described.*

The following points *interfere* with the *measurement results* and their output and should therefore *only be carried out by trained specialist personnel*. *If you are unsure about the configuration you require, please contact the service department.*

#### **5.3.2 Method**

##### **5.3.2.1 Choose**

*This menu item depends on the sensor used.*

*All analytes that can be measured* are displayed as buttons. *By touching* the buttons, the *respective method* can be switched on (i.e. *the analyte* is measured) or off (i.e. *the analyte* is not measured). *The analytes are measured simultaneously for each sample.*

Please note: light button means "on", dark button means "off".

For glucose measurement you can also select whether the *measured values* should be given as full blood values (button "*GL-Plasma-rel off*") or as plasma values (button "*GL-Plasma-rel on*" Button).

### 5.3.2.2 Adjust

The following settings can be made *here*:

- *Sample-Rerun on / off, works together with the item "sample warning limits", i.e. a sample repetition takes place if the measured value is outside the set limits*
- *Sequence number (1-999) / Daily number (starts daily at 1), determines the counting mode of the sample number*
- *Set Sample-Number, sets the sample number to the entered number*
- *Set Rack-Number, sets the cassette number to the number entered*
- *Rack-ID incr. / fix, determines whether the cassette is to be counted or whether the cassette number is to be retained*

### 5.3.2.3 Program

*Under this item, each method can be programmed individually:*

#### Unit of measurement

*For the units of measurement you can choose between mmol/l and mg/dl. For haemoglobin, g/dl is also available.*

#### Limits

*A lower and an upper warning limit can be set here. Samples found outside these limits are marked accordingly. In addition, an automatic repeat measurement can be programmed for these.*

#### Value-Correction

*Set measured value corrections*

*This function provides an analyte-dependent correction function with the equation  $y = mx + n$  to adjust the instrument to a larger laboratory instrument, if required.*

### 5.3.2.4 Show parameters

*Under this menu item, all device parameters and their configuration can be viewed and printed out.*

## 5.3.3 Functions

### 5.3.3.1 Adjust

#### EDP

*Enables the configuration of the EDP interface. This should only be done with knowledge of the EDP instruction manual.*

#### Printer TPU-80 / DPU-414

*The choice of printer must correspond to the printer connected to SUPER GL compact. In case of non-compliance, the signals sent by SUPER GL compact cannot be processed by the printer!*

#### Time

*Adjust time*

### **Language**

Currently, the dialogue languages German, English, Russian, Czech, *Polish and Slovak can be selected*. Other languages are possible on request.

### **Autocal off/on**

“Autocal on” means that after a set time regime calibrations or checks on calibrations are carried out independently of whether samples are to be measured or not. This regime is sensible if urgent samples are measured and time is critical.

“Autocal off” means that after the above time regime, calibrations lapse. In this case, a calibration is then repeated only if a sample or sample series is to be measured. If no measurements are performed for more than 6 hours, a calibration must be displayed by positioning a calibration cup. Make sure that calibration cups which have started *must not be used longer than 12 hours*.

### **Operator ID**

*The operator ID is an individual operator number for the unit user and can be activated (ID active) or deactivated (ID inactive). Touching the area changes the respective button.*

*The operator number may have a maximum of 13 digits and can be entered alphanumerically.*

*If the operator ID is activated, there are further setting options:*

- ID control / ID sample (query time of the ID: either only during control measurement or before each sample or series)
- ID scanner / ID keyboard (if a barcode scanner is connected, the *input* is requested directly via the scanner, otherwise the *input is made* via the touch keyboard)

### **Format date**

Allows entry of different date formats, e.g. European, American, etc.

*A changing of the date is not possible. If necessary, please contact our service hotline.*

### **Year**

Allows the year *to be selected* as a two or four *digit number*

### **hh:mm**

*Output format of the time, allows optional output in 24h mode or in 12h mode (AM; PM)*

### **Date-Time**

*Allows you to select the order in which the date and time are output*

### **Summer/Winter**

Allows switching on or off *for* automatic summer or winter time setting.

### 5.3.3.2 Actions

#### Washing

Rinsing of the system

*By the help of a pump solution will be pumped through the system from the suction hose in a defined time. This can be system solution or another solution, e.g. for disinfection. In this case the sensor will be destroyed.*

#### Emptying

*By the help of the pump the system will be emptied. Follow the instructions on the display.*

#### Change Syringe

*Used for maintenance (chapter 6.2)*

*Moves the plunger in a position for changing*

#### Change Needle

*Used for maintenance (chapter 6.2)*

*Moves the sample canula to a position where it is possible to change hoses and the washing chamber*

#### Show Sensor

Indicates the data currently programmed in the sensor. The display shows the following data:

Method	Method display GLU/LAC/HB)
Sensor No.	XXXXXXXXXXXX
Operating time	X months
Expiry Month	month specification
Expiry Year	year specification
Days remaining	XXX

#### Show Supply

The display shows the following data:

*OperationCard ok*

<i>cups</i>	<i>XXXX (cups mean reaction cups)</i>
<i>Month for Maintenance</i>	<i>month specification</i>
<i>Year for Maintenance</i>	<i>year specification</i>

### 5.3.3.3 Service

*Only for trained service personnel. Description in the service documentation.*

## **5.4. Control management**

### **5.4.1 General**

*To ensure effective quality assurance, SUPER GL compact also has a separate quality control function. This enables a clear tabular and graphical evaluation of the quality control process and thus the documented proof for given legal requirements. In addition, the SUPER GL compact has the option of locking the device itself for sample measurement if no control sample measurement has been carried out or if this is faulty.*

*To be able to use this quality control function, it must be activated according to chapter 5.4.2.2. Only after activation, the measured values are kept in a separate control memory in addition to the measured value memory. Only one value per user ID can be recorded daily in the control memory. The control memory can be viewed by method and position and can be statistically evaluated. The control memory holds a maximum of 100 values per control item.*

*Chapter 4.6.4 already described which setting options you have for the respective control position.*

### **5.4.2. Control management in accordance with RiliBÄK**

*RiliBÄK is the Federal Medical Council guideline on quality assurance of medical laboratory tests.*

*The quality control function is suitable for demonstrating compliance with the requirements applicable in Germany according to RiliBÄK.*

*Care must be taken that the control limits for programming the control item are identical to the method-specific limits specified in the RiliBÄK.*

*The current version of the RiliBÄK can be viewed or downloaded at [www.bundesaerztekammer.de](http://www.bundesaerztekammer.de).*

### **5.4.3 Programming of the quality control functions**

This point is reached via:

Standby -> Menue -> Controls -> Programm QC

To access this menu point you need to enter a code. This code retrieval ensures against unauthorised or unintended use. *The code is 120.*

*After entering the code, the other menu items will appear, which are described below.*

*For a better understanding, you also need the explanation of the validity of a control measurement. A control measurement is valid if the measured control value is within the control limits stored in the device and is not older than 16 hours.*

### Activate method

The control *memory* function for each of the three methods can be activated or deactivated here. This is done by pressing the relevant fields on the touch screen. *Buttons with a light background are activated.*

### Check

*Under this item, the number of measured values already in the memory is displayed for the previously activated methods, assigned to 1 and 2. By pressing the relevant field the values for the selected controls are deleted **immediately**.*

### Print protocol

*The protocol refers exclusively to the values in the control memory that have been manually changed in their validity. They are printed out in the form of a list and serve to document the changes. You can recognise these values in the display of the unit by the purple-coloured numbering.*

### Delete

*The complete control memory is deleted.*

### Save first/last Control

Here you can choose whether the first or last control measured on the day in question is recorded in the statistics (retrospective evaluation of a control cycle).

### Device-lock off/on

*If it is required that samples are not measured until a valid control measurement value is available, the device lock is switched on here.*

*The device automatically blocks all sample measurements*

- *as long as no valid control value is measured*
- *after expiry of the time-limited validity of the control measurement and*
- *when an invalid control measurement occurs*

*This function is independent of the activation of the control memory function.*

### Unlock device

*If the device lock is activated and there is no valid control measurement, this function allows you to disable the device lock for a maximum period of 16 hours. If there is another invalid control measurement, the unit is locked again.*

***The use of this function could lead to undetectable faulty sample measurements. In the case of multiple users, all affected users must be informed immediately in order to urgently evaluate the sample measurements for plausibility themselves.***

#### **5.4.4 Interpretation of the control memory**

*The control evaluation is only possible via the control memory. It is only available if the corresponding control memory function for at least one method has been activated beforehand according to 5.4.3.*

*The menu item Controls contains the sub-item Memory of controls. If you select this, you can view the control measurements for each method individually.*

*After selecting the method, all stored control values are available in the display in tabular form. By pressing the " Chart" button, the display of the control measurements can also be shown as a diagram. The display always refers to only one control position. To change the display to the second control position, touch the button " Co1/Co2".*

*In addition to the diagram, the graphical representation in the display contains the following statistically calculated values:*

- mean value of the control measurements*
- relative quadratic mean of the measurement deviation (rel. QMM)*
- coefficient of variation (CV)*

*Only valid control measurement values are used for the calculation.*

*The validity of the control measurement values can be changed manually. Basically, all measured values in the tabular display are valid. Touching the individual measured value in the table marks it as invalid. The value is indicated in the display by red digits. By touching it again, the measured value can be made valid again and gets the black digits back. All manually changed values can be identified in the display by the purple numbering.*

*When viewing the control measurements, a printer symbol always appears in the display. Touching the corresponding field on the touch screen initiates the printing process. Regardless of what is shown in the display, the statistically calculated values are always output to the printer. Invalidly marked values are then marked with "!".*

*After evaluating all control values, we recommend deleting the control memory according to 5.4.3.*

## 6 Maintenance and troubleshooting

### 6.1 Introduction

This chapter *gives information* on maintenance of the SUPER GL compact and about *problems that* may occur and how you *might be able to solve them* yourself.

**If you are unsure about *certain* aspects, DO NOT try anything you might think appropriate without qualified technical help. DO NOT open the device without an authorised service technician! Please contact our service hotline!**

### 6.2 Maintenance

The SUPER GL compact needs *maintenance* once a year by *manufacturer-authorized and trained personnel*. *Authorised personnel are people who have gained expert knowledge by completing training courses offered by the manufacturer or authorised companies.*

The touch screen *will display a message after expiration of the maintenance rate.*

Please contact the manufacturer or *distributor immediately to make an appointment for service.*

### 6.3 Servicing

The following operations can and should be carried out by the *operator*.

These *actions are part of diligent care* and serve to enhance the device's life span. *They are NOT maintenance or service work, these may only be carried out by authorised service personnel!*

#### 6.3.1 Cleaning and disinfection

Please *adhere* to the regulations valid in your laboratory with regard to cleaning and disinfecting *the device*. For disinfection wipe the entire accessible surface of the *device* with a cloth *containing* disinfectant. Use a disinfectant for surface disinfection! *Also note the instructions of the manufacturer of the disinfectant.*

## 6.3.2 Exchange of the sensor

### Removal of the sensor

Open the sensor *holder* by pulling the lever


Fig. 6.1 Sensor holder closed

Remove the sensor


Fig. 6.2 Sensor holder open

### Installation of sensor

1. Open the packaging and *take out* the sensor
2. Open the sensor *holder* by pulling the lever (fig. 6.1)
3. Insert sensor
4. Close the sensor *holder*

## 6.3.3 Exchange of the Containerkit

*The replacement of the container kit is described in chapter 3.4.1 "Containerkit SUPER GL".*

To avoid disruption in the work process of the analyser and soiling, we recommend changing the Containerkit only when the operating status is "Standby" and when there is an inbuilt *washing chamber* and sample canula. Changing the Containerkit should be carried out quickly, as liquid is expelled during the rinse processes carried out from time to time.

### **6.3.4 Decommissioning of the device / Disposal**

To decommission the device for a longer period (*more than five days*) or for transport, proceed as follows:

1. Rinse the *device* with distilled water.
2. Empty the system by *removing the hose from the supply bottle and repeating the rinsing process*.
3. Switch off the *device and disconnect*.
4. *Please store all consumables as written at the packages.*

Disposing of the device:

For the disposal of *devices* please ask your *distributor*.

### **6.3.5 Prepare *device* for transport**

To prepare the analyser for transport, proceed as follows:

- Rinse and empty system
- Switch off *the device*
- Remove Containerkit
- Remove all cups from the analyser
- Remove plug from power supply *adapter* and printer and/or EDP
- Disinfect outside of analyser!
- Pack and transport the *analyser* in the original cardboard box with foam insert, as this is the only way to avoid damage during transport as far as possible
  - Place the device in one half of the foam, making sure that the analyser is placed correctly in the mould
  - Place the other half of the foam on top and close the carton
  - *Attention, please!* If force is needed to close the transport crate, *the device* is not correctly seated in the mould.

***A non conform preparation for transport can result in expensive damages at the device.***

Please contact your *distributor* or the manufacturer if you have any further questions.

## 6.4 Error messages / Troubleshooting

### 6.4.1 Warnings

Before *printing the measured results*, the device checks whether set warning limits have been exceeded. The sample warning limits (see chapter 5.3.2.3) *are relevant for the determination of samples*, control limits (chapter 4.6.4) *are relevant for the determination of control samples*. The following warning are displayed and printed:

Warning	Explanation
++++	Values above measuring range limit
-----	Values below measuring range limit
!!	Below or above the sample warning or control limits
*!	Previous control measurement outside the control limits and below or above the sample warning

Fig. 6.3 Table of displayed and printed warnings

### 6.4.2 Errors while measuring

Zeroline unstable	<ul style="list-style-type: none"> <li>- appears incidentally -&gt; repeat incident</li> <li>- changing reagent -&gt; choose "Washing"</li> <li>- system obstructed -&gt; rinsing by hand</li> <li>- sensor defective -&gt; change sensor</li> <li>- electronic error -&gt; call service</li> </ul>
Calibration error value too low	<ul style="list-style-type: none"> <li>- choosing wrong cup</li> <li>- sensor defective -&gt; change sensor</li> <li>- electronic error -&gt; call service</li> </ul>
Calibration error drift too high	<ul style="list-style-type: none"> <li>- appears incidentally -&gt; calibrate new</li> <li>- rapid change of temperature (air draft) -&gt; calibrate new</li> <li>- system obstructed -&gt; rinsing by hand</li> <li>- sensor defective -&gt; change sensor</li> <li>- electronic error -&gt; call service</li> </ul>
Cup empty!	Remaining quantity in the reaction cup is not enough for a measurement
Sidemax left / right	<ul style="list-style-type: none"> <li>- a disturbance of the measurement signal was detected by the device</li> <li>- appears with individual samples -&gt; the sample contains substances that interfere with the measurement -&gt; measurement of this sample not possible</li> <li>- frequent incidence with samples or during calibration -&gt; sensor defective -&gt; replace sensor</li> </ul>

Fig. 6.4 Table of errors during the measuring

### 6.4.3 Device errors

Device errors caused by malfunctions of the device hardware are shown as follows:

ERROR (name of the assembly)  
(Fault description)

If such an error occurs, the user cannot act without the help of the service. The transmission of the exact error description is very important for the service for the error assessment.

## 7 Technical data

<b>Performance Data</b>	
Measuring Time	
Single Measurement Glucose / Hb	ca. 30 seconds
Single Measurement Glucose / Lactate / Hb	ca. 45 seconds
Series Measurement Glucose / Hb	ca. 30 seconds - 120 samples/hour
Series Measurement Glucose / Lactate/ Hb	ca. 45 seconds - 80 samples/hour
Measuring Range	
Glucose	0.5 mmol/l to 50 mmol/l or 9.01 mg/dl to 900.8 mg/dl
Lactate	0.5 mmol/l to 30 mmol/l or 4.5 mg/dl to 270 mg/dl
Hemoglobin	1.2 mmol/l bis 20 mmol/l bzw. 1.9 g/dl bis 32 g/dl
Precision (18 samples)	<i>Variation Coefficient Indication</i>
Glucose (12 mmol/l or 216 mg/dl)	< 1.5 %
Lactate (10 mmol/l or 90 mg/dl)	< 1.5 %
Hemoglobin (7.6 mmol/l or 12 g/dl)	< 1.5 %
Sample Volume	10 µl sample diluted with 500 µl hemolysate solution
<b>Technical Specifications</b>	
Working Temperature	+ 15 °C bis + 30 °C
Photometer Wavelength	530 nm
Supply Voltage	12 V DC
Power Consumption	max. 12 W
Classification according to Regulation 2017/746 on IVD medical devices (IVDR)	<i>In vitro diagnostic medical device</i> Risk class A
Interfaces	
Printer	serial, RS 232
EDP	serial, RS 232
Barcode reader	external, PS2
Dimensions / Weight	
Width	200 mm
Hight	170 mm
Depth	185 mm
Weight	1.7 kg

Fig. 7.1 List of technical data

*\*The performance data is largely determined by the sensors used and the other consumables. The data given here represents the technical capabilities of the device under optimal conditions.*