
# HemoCue® Glucose 201 **RT**

Operating Manual

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**EN****Table of Contents**

Introduction.....	3
Components.....	4
Start-up.....	6
Set-up: Audio signal, time and date.....	8
Measuring: Capillary blood.....	10
Measuring: Control material, venous or arterial blood.....	18
Set-up: QC Test.....	20
Set-up: Memory function – Scroll.....	22
Set-up: Delete results.....	24
Set-up: Printer and connectivity.....	26
Maintenance: Daily maintenance.....	28
Maintenance: Optical parts.....	30
Troubleshooting Guide.....	32
Specifications.....	48

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# HemoCue® Glucose 201 RT System

EN

Thank you for choosing the HemoCue® Glucose 201 RT system. The system is a fast, simple and reliable method for the determination of glucose in whole blood. The results obtained are comparable to those from standard laboratory methods. Capillary, venous or arterial whole blood may be used. The system is intended for use in near-patient (point-of-care) and laboratory testing. The system is for *In Vitro* Diagnostics use only. The HemoCue Glucose 201 RT Analyzer is available with or without plasma conversion and only to be used with HemoCue Glucose 201 RT Microcuvettes. Analyzers with plasma conversion will show P-equiv on the display.

*Read and follow this operating manual and the HemoCue Glucose 201 RT Microcuvettes package insert to attain optimum performance and safety. Any other use of the system than recommended by the manufacturer may impair the safety.*

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
Components  
Bestandteile  
Composants  
Componenten

1. HemoCue Glucose 201 RT Analyzer \*
2. Power adapter \*\*
3. 4 type AA, 1.5 V Alkaline batteries
4. HemoCue Glucose 201 RT Microcuvettes
5. HemoCue Glucose 201 RT Operating Manual and HemoCue Glucose 201 RT Quick Reference Guide

#### Materials provided

HemoCue Glucose 201 RT Analyzer  
Power adapter  
HemoCue Glucose 201 RT Operating Manual  
HemoCue Glucose 201 RT Quick Reference Guide

The HemoCue® Glucose 201 RT Analyzer and its power adapter are delivered in a carton. Open the carton on a stable surface and lift out the analyzer and other components. If no AC power is available, use 4 type AA Alkaline batteries, 1.5 V. On the bottom of the analyzer there is a lid covering the battery compartment. Press the flap to remove the lid. Place the batteries in the battery compartment and replace the lid. Consult local environmental authorities for proper disposal of batteries.

\* Do not open the cover of the analyzer. Note: The warranty is voided if the analyzer has been opened.

\*\*  Only use adapters as listed under Specifications.

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
- 1a. Attach the power adapter to the power inlet at the back of the analyzer.
- 1b. If no AC power is available, insert four type AA Alkaline batteries, 1.5 V, into the battery compartment.
- 1c. If a battery symbol appears on the display, the batteries are running low on power. The analyzer will continue to give accurate results but the batteries should be replaced as soon as possible.
2. The analyzer can be connected directly to a printer (see section *Set-up Printer and connectivity*) or a PC (contact your local distributor or HemoCue AB for further information).
3. Pull the cuvette holder out to the loading position.
4. Press and hold the left button (4a) until the display is activated (all symbols appear on the display) (4b)\*. The analyzer performs an internal quality control, the self-test and after approximately 10 seconds the display will show three flashing dashes and the HemoCue symbol (4c). This indicates that the analyzer is ready for use.

To turn the analyzer off, press and hold the left button until the display reads OFF and then goes blank. Disconnect the analyzer from the power source by removing the power adapter from the outlet.

\* An analyzer without plasma conversion will not show "P-equiv" on the display.

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






1. Press both buttons at the same time.
2. The display now shows a flashing QC-symbol.
3. Use the right button to scroll until the display shows an audio symbol in the upper right corner (3a). The signal can be activated (3a) or deactivated (3b) by pressing the left button.
4. Continue to scroll using the right button until the display shows characters for time, date and year. The hour figure will be flashing.
5. To change the hour, either hold the left button down for faster advancement or press and release the left button to advance slowly. When the desired hour is displayed, press and release the right button. The minute figure will then be flashing. Proceed as above until desired time, year and date are obtained.
6. When all the settings are completed, hold the right button down for approximately 5 seconds. The analyzer will then automatically return to the measuring position.

The analyzer is supplied with default time and date settings.

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
**Mesure** Sang capillaire

**Meting** Capillair bloed

 Always wear protective gloves. Handle blood with care, as it may be infectious. Follow local safety procedures for disposal of used microcuvettes. Operating temperature: 15–27 °C (59–80 °F). Start measurement no later than 40 seconds after filling the microcuvette.

Note: In cases of severe hypotension or peripheral circulatory failure, glucose measurement from capillary samples may be misleading. In such circumstances it is recommended that the glucose level should be measured using venous or arterial whole blood.

1. To perform a test using capillary blood, the cuvette holder should be in its loading position. The display will show three flashing dashes and the HemoCue symbol.
2. Make sure the patient's hand is warm and relaxed. Use only the middle or ring finger for sampling. Avoid fingers with rings on.
3. Clean the finger with disinfectant and allow to dry completely or wipe off with a dry, lint-free wipe.
4. Using your thumb, lightly press the finger from the top of the knuckle towards the tip.

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
13. Place the microcuvette into the cuvette holder and start measurement as soon as possible but no later than 40 seconds after filling the microcuvette, by gently sliding the cuvette holder to the measuring position.
14. During the measurement "Ⓜ" and three fixed dashes will be shown.
15. Most results will appear on the display within 1 minute for glucose levels  $\leq 7$  mmol/L ( $\leq 126$  mg/dL). The result will remain on the display as long as the cuvette holder is in the measuring position. An analyzer with plasma conversion will show "P-equiv" on the display. If not in use, an analyzer operated on battery power will automatically turn off after 5 minutes. Do not remeasure the cuvette.
16. Discard the microcuvette after measurement. Microcuvettes are for single use only. Follow local safety procedures for disposal of used microcuvettes.

All the results are stored in the analyzer and can be viewed on the display using the scroll function (see section *Set-up Memory function-scroll*) or by using a PC program (contact your local distributor or HemoCue AB for further information).

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
**Measuring** Control material, venous or arterial blood

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**Meting** Controlematerial, veneus of arterieel bloed

**!** Always wear protective gloves. Handle blood with care, as it may be infectious. Follow local safety procedures for disposal of used microcuvettes.  
Operating temperature: 15–27 °C (59–80 °F). Start measurement no later than 40 seconds after filling the microcuvette.

1. If the specimen has been stored in a refrigerator, allow it to warm up to room temperature before mixing. The blood or control material should be well mixed prior to measurement.
2. Place a drop of blood or control material (>4 µL) onto a hydrophobic surface using a pipette (2a) or other suitable transfer device (2b).
3. a. Filling end; b. Optical eye; c. Air hole.
4. Fill the microcuvette at the filling end (3a). It is important to fill the microcuvette in one continuous process. Do not refill! If a second sample is to be taken, fill a new cuvette from a new drop of specimen. This should not be done until the measurement of the first sample is completed.
5. Wipe off excess specimen from the outside of the microcuvette with a clean, lint-free wipe being careful not to touch the filling end (3a) of the microcuvette, which could result in specimen being drawn out of the microcuvette.
6. Visually inspect the microcuvette. If the microcuvette is not completely filled with specimen or if there are air bubbles, discard and fill a new microcuvette. Small bubbles around the filling end (3a) can be ignored.

Perform the analysis as per steps 13–16 in section *Measuring Capillary blood*.

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
### Configuration test CQ

## Set-up QC-test

To specifically identify QC results in the memory, proceed as follows:

1. Press both buttons at the same time.
2. The display shows a flashing QC-symbol.
3. Select QC-test by pressing the left button.
4. The analyzer automatically returns to its measuring mode and the QC- symbol appears on the display.
5. Fill the cuvette with the control solution and carry out the measurement in accordance with section *Measuring Control material, venous and arterial blood*.
6. Once the QC measurement has been made, the analyzer automatically returns to the measuring mode and the QC-symbol disappears from the display. Patient sampling can begin.

If more QC-tests are to be performed, proceed in the same way following steps 1–6. To deactivate the QC-test, follow steps 1–2 and scroll using the right button until another set-up activity is shown on the display. Hold the right button down for approximately 5 seconds. The analyzer automatically returns to its measuring position and the QC-symbol is no longer visible.

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
Set-up Geheugenfunctie – Scrollen

The analyzer automatically stores up to 600 results. When the memory is full the analyzer will automatically overwrite the oldest result. To view results stored in the memory use the following procedure.

1. The cuvette holder may be in the loading position (1a) or in the measuring position (1b). If in the measuring position the most recent measurement is displayed (1b).
2. Press the left or right button to scroll backwards or forwards between the results. The stored values are now shown on the display. When the button is held down, the user can see the order of the results (the most recent measurement is always number 1) and when the button is released the result is shown.
3. If no buttons are pressed while the cuvette holder is in the loading position (3a), the analyzer automatically resets after 5 seconds and new measurements can be made. If in the measuring position (3b), pull the cuvette holder out to the loading position (3a) and wait for the HemoCue symbol. The analyzer is then ready for a new measurement.

Note: The analyzer will not respond to other commands while a measurement is being made. If the error code E35 is shown while scrolling through the results, the measurement has not been stored in the memory, see the *Troubleshooting Guide*.

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
1. Press both buttons at the same time.
2. The display now shows a flashing QC-symbol.
3. Scroll using the right button until the display shows a flashing waste bin in the left corner.
4. To delete the current result, press and hold the left button. A warning triangle together with a "I" appears on the display (4a) for approximately 3 seconds. The current result has now been deleted and the display shows a "0" (4b) to confirm that this has been done. Only the current result can be deleted, which must be done immediately after the measurement has been made.
5. To delete all stored results, hold both buttons down at the same time for 5 seconds. The display shows the waste bin, the warning triangle and the total number of results stored (5a). Continue to hold both buttons down until the number changes to "000" and the warning triangle disappears (5b). All stored results have now been deleted.

Note: Deleted results cannot be recovered. If results are to be saved, they must be printed out or downloaded to a PC before starting the deleting process.

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
Set-up Printer en connectiviteit

Only use ASCII printers recommended by HemoCue. Always follow instructions for use for the printer.

1. Connect analyzer to printer or computer\* via serial or USB port (analyzer can be turned on or off).
2. When measurement is completed, the result will automatically be printed or transferred to the computer with date and time.

Time and date must be activated before data can be printed.

\* The analyzer can be connected to HemoCue software or external software applications for transfer of data, using either a USB cable or serial cable. For more information, contact your local distributor or HemoCue AB.

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

Maintenance Daily maintenance

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Maintenance Maintenance quotidienne

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The cuvette holder should be cleaned after each day of use.

1. Turn off the analyzer. Pull the cuvette holder out to the loading position.
2. While pressing the catch, carefully rotate the cuvette holder sideways in open position as far as possible to the left.
3. Remove the cuvette holder from the analyzer, it will come off the stainless steel pin it rotates on.
4. Clean the cuvette holder with alcohol (20–70%) or mild detergent

Wait 15 minutes before replacing the cuvette holder and using the analyzer. Make sure the cuvette holder is dry before inserting.

The optical parts should be cleaned when directed to do so in the *Troubleshooting Guide*, see section *Maintenance Optical parts*.

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



Maintenance Optical parts

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Dirty optical parts may cause an error code. Follow step 1–3 under *Daily maintenance* section and then clean optical parts as follows.

Note! Make sure that the HemoCue® Cleaner reaches both upper and lower cover glasses, see pictures 5a and 5c.

5. With the cuvette holder removed from the analyzer push the Cleaner into the opening of the optic unit, as far in as possible. Move from side to side 5–10 times (5b) and thereafter, placing the Cleaner to the left, push in and pull out the Cleaner 5–10 times (5c) to clean the cover glasses. If the Cleaner is stained, repeat with a new Cleaner.
6. Wait 15 minutes before replacing the cuvette holder and using the analyzer. Make sure the cuvette holder is dry before inserting it.

The cover may be cleaned with alcohol (20–70%) or mild detergent. Note! As an alternative to the HemoCue® Cleaner, a lint free cotton swab, non-pretreated, moistened with alcohol (20–70% without additive) or water may also be used. If a cotton swab is used make sure it is not too wet and nor too dry. Use a dry swab to wipe away excess liquid in the optic house after cleaning with a moistened swab. To avoid scratches on the cover glasses, only the cotton part of the swab should come in contact with the cover glasses.

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# Troubleshooting Guide

If you are unable to resolve a problem by following this Troubleshooting Guide, please contact your local distributor or HemoCue AB. The analyzer should be cleaned prior to service, see section *Maintenance* for instructions. The analyzer has no serviceable parts. Note! Do not open the cover of the analyzer. The warranty is voided if the analyzer has been opened.

Symptom	Explanation	Action
The analyzer shows an error code	May be a temporary fault.	Turn off the analyzer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below.
E00	No stable endpoint of the measurement is found within the time range. 1. The cuvette is faulty. 2. The circuit board is out of order.	1a. Check the expiration date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The analyzer needs service. Contact your local distributor.
E01–E05	1. Dirty optronic unit or faulty electronic or optronic unit.	1a. Turn off the analyzer and clean the optronic unit as described in the <i>Maintenance</i> section. 1b. The analyzer needs service. Contact your local distributor.
E06	1. Unstable blank value. The analyzer might be cold.	1. Turn off the analyzer and allow it to reach operating temperature. If the problem continues, the analyzer needs service. Contact your local distributor.
E07	1. The battery power is too low.	1a. The batteries need to be replaced. Turn off the analyzer and replace the batteries (4 type AA, 1.5 V). 1b. Use the power adapter.
E08	The absorbance is too high. 1. Light blocking item in the cuvette holder.	1a. Check that the analyzer and microcuvettes are being used according to the HemoCue Glucose 201 RT operating manual and the package insert. 1b. The analyzer needs service. Contact your local distributor.
E09–E30	1. Dirty optronic unit or faulty electronic or optronic unit.	1a. Turn off the analyzer and clean the optronic unit as described in the <i>Maintenance</i> section. 1b. The analyzer needs service. Contact your local distributor.

Symptom	Explanation	Action
E35	<ol style="list-style-type: none"> <li>The current measurement has been incorrectly stored and will remain as error code E35 until the entire memory has been deleted.</li> </ol>	<ol style="list-style-type: none"> <li>No action is normally required.</li> <li>If storage problems continue, the analyzer needs service. Contact your local distributor.</li> </ol>
E70–E71	<p>Erroneous kinetics.</p> <ol style="list-style-type: none"> <li>The cuvette is faulty.</li> <li>The sample might be grossly lipemic.</li> </ol>	<ol style="list-style-type: none"> <li>Check that the system is used according to the instruction for use. Fill a NEW cuvette and perform a measurement.</li> <li>If the error code appears again, use a suitable laboratory reference method to analyze the specimen.</li> </ol>
ECU	<ol style="list-style-type: none"> <li>The cuvette is empty or not filled with enough sample. Please note, the ECU function captures primarily empty cuvettes, not under-filled cuvettes.</li> <li>No chemical reaction is identified in the blood filled cuvette.</li> </ol>	<ol style="list-style-type: none"> <li>Fill a NEW cuvette and perform a measurement.</li> <li>If the ECU function appears again, contact your local distributor or HemoCue AB. Use another HemoCue analyzer or a suitable laboratory reference method to analyze the specimen.</li> </ol>
HHH	<ol style="list-style-type: none"> <li>Measured value exceeds 27.8 mmol/L (500 mg/dL) for a whole blood system and 30.8 mmol/L (555 mg/dL) for a plasma equivalent system.</li> </ol>	<ol style="list-style-type: none"> <li>The result is above the measuring range.</li> </ol>

Symptom	Explanation	Action
No characters on the display	<ol style="list-style-type: none"> <li>1. The analyzer is not receiving power.</li> <li>2. If on battery power, the batteries need to be replaced.</li> <li>3. The display is out of order.</li> </ol>	<ol style="list-style-type: none"> <li>1a. Check that the power adapter is connected to the AC power supply.</li> <li>1b. Check that the power adapter is securely connected to the analyzer.</li> <li>1c. Check that the cable is not damaged.</li> <li>2. Turn off the analyzer and replace the batteries (4 type AA, 1.5 V).</li> <li>3. The analyzer needs service. Contact your local distributor.</li> </ol>
The display gives erroneous characters	<ol style="list-style-type: none"> <li>1. The display is out of order.</li> <li>2. The microprocessor is out of order.</li> </ol>	<ol style="list-style-type: none"> <li>1. The analyzer needs service. Contact your local distributor.</li> <li>2. The analyzer needs service. Contact your local distributor.</li> </ol>
The display shows 	<ol style="list-style-type: none"> <li>1. The batteries need to be replaced.</li> <li>2. If on AC power, the power adapter or the circuit board is out of order.</li> </ol>	<ol style="list-style-type: none"> <li>1. Turn off the analyzer and replace the batteries (4 type AA, 1.5 V).</li> <li>2a. Check that the power adapter is properly connected and working.</li> <li>2b. The analyzer needs service. Contact your local distributor.</li> </ol>
The display does not switch from  and "GLU" to three flashing dashes and  (ready for measuring)	<ol style="list-style-type: none"> <li>1. The magnet in the cuvette holder may be missing.</li> <li>2. The magnetic sensor is out of order.</li> </ol>	<ol style="list-style-type: none"> <li>1. The analyzer needs service. Contact your local distributor.</li> <li>2. The analyzer needs service. Contact your local distributor.</li> </ol>

Symptom	Explanation	Action
Too high or low values for controls or patient samples compared to expected values	<ol style="list-style-type: none"> <li>1. The microcuvettes are damaged, beyond their expiration date or have been improperly stored.</li> <li>2. The optical eye of the microcuvette is contaminated.</li> <li>3. The control has not been properly mixed and/or has not reached room temperature.</li> <li>4. Air bubbles are present in the cuvette.</li> <li>5. The optronic unit is contaminated.</li> <li>6a. The control used is not suitable for the HemoCue Glucose 201 RT system.</li> <li>6b. The control is beyond its expiration date or has been stored incorrectly.</li> <li>7. Incompletely filled microcuvette.</li> <li>8. The measurement needs to be started no later than 40 seconds after filling the microcuvette.</li> <li>9. Improper sampling technique.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check the expiration date and the storage conditions of the microcuvettes.</li> <li>2. Measure the control/sample with a new cuvette.</li> <li>3. Make sure that the control is properly mixed and that it has reached room temperature. If the problem continues, contact the manufacturer of the control.</li> <li>4. Check the microcuvette for air bubbles. Remeasure the control/sample with a new microcuvette.</li> <li>5. Clean the optronic unit as instructed under <i>Maintenance</i> section.</li> <li>6a. If a quality control test is to be performed, only use quality controls recommended by HemoCue, see relevant package insert for more information.</li> <li>6b. Check the expiration date and the storage conditions of the control. Repeat the measurement with a new control/sample. If the problem continues, contact the manufacturer of the control.</li> <li>7. Measure the control/sample with a new microcuvette.</li> <li>8. Measure the control/sample with a new microcuvette.</li> <li>9. Follow instructions in section <i>Measuring</i>.</li> </ol>
Measurements on patient samples are higher or lower than anticipated	<ol style="list-style-type: none"> <li>1. Improper sampling technique.</li> <li>2. The cuvettes are beyond their expiration date, damaged or have been improperly stored.</li> <li>3. The optical eye of the cuvette is contaminated.</li> <li>4. Air bubbles in the cuvette</li> <li>5. The optronic unit is dirty.</li> <li>6. The calibration of the analyzer has changed.</li> </ol>	<ol style="list-style-type: none"> <li>1. Follow instructions in section <i>Measuring</i>.</li> <li>2. Check the expiration date and the storage conditions of the cuvettes.</li> <li>3. Remeasure the sample with a new cuvette.</li> <li>4. Check the cuvette for air bubbles. Remeasure the sample with a new cuvette.</li> <li>5. Clean the optronic unit as described in the <i>Maintenance</i> section.</li> <li>6. The analyzer needs service. Contact your local distributor.</li> </ol>

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# Specifications

## Intended purpose

The HemoCue Glucose 201 RT System is intended for quantitative determination of glucose in whole blood (capillary, venous or arterial), supplementing the clinical evidence in the diagnosis of diabetes, the detection of idiopathic hypoglycemia, monitoring patients with diabetes and monitoring neonatal blood glucose levels. The system should not be used on critically ill neonates in neonatal intensive care settings. The HemoCue Glucose 201 RT System is an automated system for professional use by healthcare professionals and laboratory personnel, intended for near-patient (point-of-care) and laboratory testing. The HemoCue Glucose 201 RT System is for *In Vitro* Diagnostics use only. The HemoCue Glucose 201 RT Analyzer is only to be used with HemoCue Glucose 201 RT Microcuvettes.

## IVD Medical Device Regulation

The HemoCue Glucose 201 RT System complies with the IVD Medical Device Regulation (EU) 2017/746 and carries the CE mark.

## Principles of the method/procedure

### *Principle of the method*

The chemical reaction in the cavity has two phases, hemolysis and the glucose reaction. The glucose reaction is a modified glucose dehydrogenase method in which a tetrazolium salt (MTT) is used to obtain a quantification of glucose in visible light.  $\alpha$ -D-glucose is transformed to  $\beta$ -D-glucose using mutarotase. Glucose dehydrogenase acts as a catalyst for the oxidation of  $\beta$ -D-glucose, to form NADH, which in the presence of diaphorase produces a colored formazan with MTT.

### *Principle of the procedure*

The system consists of an analyzer together with microcuvettes. The microcuvette serves both as a pipette and as a measuring cuvette and is for single-use only. A blood sample of approximately 4  $\mu$ L is drawn into the cavity by capillary action. The measurement takes place in the analyzer in which the transmittance is measured and the absorbance and glucose level is calculated. The calibration of the HemoCue Glucose 201 RT System is traceable to the ID GC-MS reference method which is calibrated by NIST Standard reference material. The HemoCue Glucose 201 RT is factory calibrated and needs no further calibration.

## Composition

The microcuvette is made of copolyester plastic and contains the following reagent: <75  $\mu$ g/g cuvette MTT (methylthiazolyldiphenyl tetrazolium), <130  $\mu$ g/g cuvette saponin, <40  $\mu$ g/g cuvette sodium fluoride, <525  $\mu$ g/g cuvette enzyme mix: mutarotase (bacterial), glucose dehydrogenase (bacterial), diaphorase (bacterial), NAD (nicotinamide-adenine-dinucleotide), and non-reactive ingredients.

## Warning and precaution

In accordance with "Good Laboratory Practices", intervention levels and normal values should be established by each laboratory when instituting a new method. The specimen collection should be done according to instruction for use to reduce the risk of incorrect results. In cases of severe hypotension and peripheral circulatory failure, glucose measurement from capillary samples may be misleading. In such circumstances it is recommended that the glucose level should be measured using venous or arterial whole blood. The measurement needs to be started no later than 40 seconds after filling the microcuvette. If a second sample is to be taken, fill a new microcuvette from a new drop of blood. This should not be done until the measurement of the first sample is completed. Do not remeasure the microcuvette. If the result exceeds the measuring range of the system "HHH" will be displayed. Always handle blood specimens with care as they may be infectious. Consult local environmental authorities for proper disposal.

## Storage and handling

### *HemoCue Glucose 201 RT Microcuvettes*

The microcuvettes are to be stored at 0–30 °C (32–86 °F). Unused HemoCue Glucose 201 RT Microcuvettes should be kept in the original package but once an individual package is opened, the microcuvette must be used immediately or discarded. Use the microcuvettes prior to the expiration date that is printed on the package.

### *HemoCue Glucose 201 RT Analyzer*

The analyzer can be stored and transported at 0–50 °C (32–122 °F). Operating temperature is 15–27 °C (59–80 °F). Allow the analyzer to reach operating temperature (app. 30 minutes) before use. The analyzer should not be operated at high (i.e. >90% non-condensing) humidity.

### Specimen handling and preparation

Capillary, venous or arterial blood may be used. To achieve optimal result quality the specimen collection technique is essential. Make sure, regardless of type of blood, to do collection according to instructions for use. For venous or arterial specimens, one of the following anticoagulants shall be used: EDTA, sodium heparin or lithium heparin with or without gel and one of the following glycolysis inhibitors with anticoagulants may be used: sodium fluoride/potassium oxalate, sodium fluoride/EDTA or sodium fluoride/heparin.

Glycolysis is a major concern in all glucose measurements. To minimize the effect of glycolysis, measure the blood sample as soon as possible after taking the sample. Samples of blood collected in containers with recommended anticoagulants should be analyzed within 30 minutes (when stored in 15–25 °C (59–77 °F)). Mix all specimen tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8–10 times by hand.

Note! Samples in lithium heparin gel tubes cannot be used after centrifugation.

Note! Plasma specimen cannot be used.

### Materials required but not provided

- HemoCue Glucose 201 RT Microcuvettes
- Lancet (capillary samples)
- Pipette or other transfer device (venous or arterial samples)
- Hydrophobic surface (for venous or arterial blood)
- Lint-free tissue (non-fraying)

### Quality control

The HemoCue Glucose 201 RT Analyzer has an internal quality control, the self-test. Every time the analyzer is turned on, it will automatically verify the performance of the analyzer. The analyzer will give an error code E01–E07 if the measuring electronics is unstable in any way. Upon passing the self-test, the display will show the HemoCue symbol and three flashing dashes, indicating that the analyzer is ready to perform a measurement.

Follow local guidelines regarding quality control procedures. If an external quality control is required by local or other regulation, only use controls recommended by HemoCue, see relevant package insert for more information.

### Reference values

Fasting glucose values (reference interval)\*:

Whole blood glucose, adults 3.6–5.3 mmol/L (65–95 mg/dL).

Plasma glucose, adults 4.1–5.5 mmol/L (74–99 mg/dL).

For diagnosis of diabetes mellitus, follow local recommendations or use the following value according to WHO and ADA:

Fasting whole blood glucose capillary or venous  $\geq 6.1$  mmol/L ( $\geq 110$  mg/dL).

Fasting plasma glucose capillary or venous  $\geq 7.0$  mmol/L ( $\geq 126$  mg/dL).

\* Rifai, Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics, 8th Ed., Elsevier, 2018.

### Limitations of the method/procedure

For limitations of the method/procedure, see HemoCue Glucose 201 RT Microcuvettes package insert.

## Specific performance characteristics

**Table 1.** Glucose 201 RT System ranges.

	Measuring range	Displayed range	Limit of Detection (LoD)	Linear range
Whole blood system	0.55–27.8 mmol/L (10–500 mg/dL)	0–27.8 mmol/L (0–500 mg/dL)	0.55 mmol/L (10 mg/dL)	0–27.8 mmol/L (0–500 mg/dL)
Plasma equivalent system	0.61–30.8 mmol/L (11–555 mg/dL)	0–30.8 mmol/L (0–555 mg/dL)	0.61 mmol/L (11 mg/dL)	0–30.8 mmol/L (0–555 mg/dL)

If the results exceed the measuring range of the system “HHH” will be displayed. If the patient assessment does not align with results from the HemoCue Glucose 201 RT System the test results should be verified against a reference method.

**Table 2.** Glucose 201 RT accuracy and bias, capillary and venous blood samples (results in whole blood).

Sample type	N	Min–Max mg/dL	Regression line	Correlation coefficient r	Glucose level mg/dL	Bias
Venous (vs. ID GC-MS)	104	32.6–460.5	Y = 1.06X – 3.55	0.998	0–50 51–120 121–200 201–300 301–400 401–500	-4.5 mg/dL 2.2 mg/dL 2.9% 4.6% 4.4% 4.3%
Capillary (vs. YSI 2300 STAT Plus)	102	13.1–476.0	Y = 1.00X + 2.15	0.997	0–120 121–400 401–500	1.3 mg/dL 1.3% 0.7%

Note: To convert values in mg/dL to mmol/L divide by 18.016

**Table 3.** Glucose 201 RT repeatability and within-laboratory precision (results for liquid control material).

Control level	N	Mean mg/dL (mmol/L)	Repeatability		Within-laboratory precision	
			SD <sub>pooled</sub> mg/dL (mmol/L)	CV %	SD <sub>pooled</sub> mg/dL (mmol/L)	CV %
1	800	45.9 (2.5)	1.15 (0.064)	2.5	1.28 (0.071)	2.8
2	800	101.8 (5.7)	1.36 (0.075)	1.3	1.62 (0.090)	1.6
4	800	305.7 (17.0)	3.26 (0.180)	1.1	3.72 (0.210)	1.2

Note: To convert values in mg/dL to mmol/L divide by 18.01

## Technical Specifications

Dimensions: 85 × 160 × 43 mm (3.35 × 6.3 × 1.69 inches)

Weight: 350 g (0.77 lbs) with batteries installed

Battery: 4 batteries type AA, 1.5 V Alkaline

Power adapter: CE marked. Type FW7333SM/12, FW7660M/12, FW8002.1/12

Only use power adapters recommended by HemoCue.

Analyzer electrical rating: 12 V  $\overline{=}$ , 1 W

Pollution degree: 2

Overvoltage category: II

Altitude: up to 2000 m above sea level

The instrument is made for continuous mode

## Warranty

The analyzer carries a 24-month warranty from the day of receipt. After the expiration date of the warranty, maintenance and repairs are offered at a fixed price. For technical difficulties or repair contact your local distributor or HemoCue AB. Any other use of the system than recommended by the manufacturer will void the warranty.

## Service and disposal

The Analyzer should be cleaned as recommended under section *Maintenance* prior to service or disposal. Dispose of the battery lid and cuvette holder as potentially infectious waste. Remove the batteries. Clean the Analyzer cover and optical parts according to section *Maintenance*. Disinfect the Analyzer external surfaces and optical parts according to local procedures. Wait for three weeks, then dispose of the Analyzer according to local procedures for electronic waste. For further information, contact HemoCue distributor or HemoCue AB.

## EMC and Electrical Safety

The analyzer has been tested and complies with following standards:

- IEC/EN 61010-1
- IEC/EN 61010-2-101
- IEC/EN 61326-2-6 (including applicable parts of IEC/EN 61326-1)

The analyzer has been tested for indoor use. It is recommended to evaluate the electromagnetic environment prior to use of the device.

Do not use this device in close proximity to sources of strong electromagnetic radiation (e.g unshielded intentional RF sources), as these can interfere with the proper operation.

*Note: It is the manufacturer's responsibility to provide equipment electromagnetic compatibility information to the customer or user.*

*Note: It is the user's responsibility to ensure that a compatible electromagnetic environment for the equipment can be maintained in order that the device will perform as intended.*

## Spare parts

- Power adapter\*
- Battery Lid
- Cuvette holder

*\* Not available in all countries*

## Optional items

- Printer
- Software applications
- Cables for PC and printer connection
- HemoCue Cleaners
- Hard Carrying Case
- HemoCue Safety Lancets

## References

- NIST Standard reference material, <https://www.nist.gov/srm>.
- Atkin et al, Annals of Internal Medicine, 1991, 114:12, 1020–1024.
- Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800.
- Rifai, Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics, 8th Ed., Elsevier, 2018.
- Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycaemia. Report of a WHO/IDF Consultation. Geneva: WHO 2006.
- Definition, Diagnosis and Classification of Diabetes Mellitus and its Complications. Report of a WHO Consultation. Geneva: WHO 1999.
- American Diabetes Association. Diagnosis and Classification of Diabetes Mellitus, Diabetes Care, 2006, 29: S43–S48.
- Interference Testing in Clinical Chemistry NCCLS Approved Guideline; CLSI Document EP7-A2 Vol. 25 No. 27.
- HemoCue Glucose 201 RT Microcuvettes Package Insert.

## Reporting of serious incidents

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to HemoCue AB or local distributor, and to your national authority as applicable according to local requirements.

## Symbols used

Caution

Consult instructions for use

Device for Near-patient testing

*In Vitro* Diagnostic medical device

Biological risk

Manufacturer and date of manufacture

Serial number

Catalogue number

Temperature limitation

Humidity limitation

Only valid within the European Community.  
Indicates separate collection for waste of  
electrical and electronic equipment.

DC inlet

USB port

Serial port

Class II equipment  
*on power adapter*

For indoor use only  
*on power adapter*

## Revision history

Rev. 220524 First revision to comply with IVDR

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## Literatur

- NIST Standard reference material, <https://www.nist.gov/srm>.
- Atkin et al, Annals of Internal Medicine, 1991, 114:12, 1020–1024.
- Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800.
- Rifai, Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics, 8th Ed., Elsevier, 2018.
- Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycaemia. Report of a WHO/IDF Consultation. Geneva: WHO 2006.
- Definition, Diagnosis and Classification of Diabetes Mellitus and its Complications. Report of a WHO Consultation. Geneva: WHO 1999.
- American Diabetes Association. Diagnosis and Classification of Diabetes Mellitus, Diabetes Care, 2006, 29: S43–S48.
- Interference Testing in Clinical Chemistry NCCLS Approved Guideline; CLSI Document EP7-A2 Vol. 25 No. 27.
- Packungsbeilage HemoCue Glucose 201 RT Microcuvettes.

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## Références

- NIST Standard reference material, <https://www.nist.gov/srm>.
- Atkin et al, *Annals of Internal Medicine*, 1991, 114:12, 1020–1024.
- Hannestad et al, *Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood*, *Clinical Chemistry*, 1997; 43:5, 794–800.
- Rifai, *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*, 8th Ed., Elsevier, 2018.
- *Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycaemia*. Report of a WHO/IDF Consultation. Geneva: WHO 2006.
- *Definition, Diagnosis and Classification of Diabetes Mellitus and its Complications*. Report of a WHO Consultation. Geneva: WHO 1999.
- American Diabetes Association. *Diagnosis and Classification of Diabetes Mellitus*, *Diabetes Care*, 2006, 29: S43–S48.
- *Interference Testing in Clinical Chemistry* NCCLS Approved Guideline; CLSI Document EP7-A2 Vol. 25 No. 27.
- Notice HemoCue Glucose 201 RT Microcuvettes.

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- NIST Standard reference material, <https://www.nist.gov/srm>.
- Atkin et al, Annals of Internal Medicine, 1991, 114:12, 1020–1024.
- Hannestad et al, Accuracy and precise isotope dilution mass spectrometry method for determining glucose in whole blood, Clinical Chemistry, 1997; 43:5, 794–800.
- Rifai, Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics, 8th Ed., Elsevier, 2018.
- Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycaemia. Report of a WHO/IDF Consultation. Geneva: WHO 2006.
- Definition, Diagnosis and Classification of Diabetes Mellitus and its Complications. Report of a WHO Consultation. Geneva: WHO 1999.
- American Diabetes Association. Diagnosis and Classification of Diabetes Mellitus, Diabetes Care, 2006, 29: S43–S48.
- Interference Testing in Clinical Chemistry NCCLS Approved Guideline; CLSI Document EP7-A2 Vol. 25 No. 27.
- Bijsluiters van HemoCue Glucose 201 RT Microcuvettes.

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© 2022 HemoCue AB
HemoCue AB | Kuvettgatan 1 | SE-262 71 Ängelholm | Sweden | Phone +46 77 570 02 10 | Fax +46 77 570 02 12